[www.truking.com](http://www.truking.com)

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No.888 Watertown Avenue, Jiutai Economic Development Area, Changchun, Jilin Province.

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The company reserves the right to interpret the contents of this manual without prior notice. The pictures in this manual are for reference and subject to the actual conditions.

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
TRUKING WATERTOWN PHARMACEUTICAL EQUIPMENT CO.,LTD.

**让世界制药工业插上智慧的翅膀**  
Equip global pharmaceutical manufacturing industry with intelligent wings.

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Truking Watertown Pharmaceutical Equipment Co., LTD., a wholly-owned subsidiary of Truking Technology Limited, was founded in 1993. It is the leading R&D and manufacturing base of pharmaceutical water system (preparation, storage and distribution), liquid preparation system, stainless steel pressure vessel, sewage treatment, evaporation, crystallization and heat exchanger in Asia. With more than 1300 employees and 2 billion yuan of assets, the company serves more than 1000 pharmaceutical enterprises around the world.

Truking Watertown is located in Jiutai Economic Development Zone, Changchun. It is a member of China Pharmaceutical Equipment Industry Association, China Pharmaceutical Equipment Engineering Association and China Membrane Industry Association. It has been rated as a national high-tech enterprise, Provincial Specialized and Sophisticated Enterprise, Provincial Enterprise Technology Center. It has been honored as provincial innovative technology Enterprise, Changchun Top 100 Private Enterprise, and Changchun Small Giant Technology Enterprise. Truking Watertown has Chinese special equipment Class D (Class I, Class II and Class III below 10MPa) pressure vessel design and manufacturing license, Pressure pipeline (GC2) installation license, construction mechanical and electrical installation engineering contracting Class II, environmental engineering contracting Class II, environmental engineering design special (water pollution control engineering) Class B, ASME(U) steel seal design and design manufacturing qualification, EU CE(PED, MD) certification qualification. Truking Watertown has applied for 320 patents, 180 patents of which have been authorized until April 2023, and one PCT international patent application. Products have been exported to the United States, Italy, Belgium, Greece, Turkey, Russia, South Korea, Thailand, Vietnam, Indonesia, India, Egypt, Saudi Arabia, Morocco, Peru and other 30 countries and regions, the international market share is rapidly increasing year by year.

Truking Watertown has an experienced R&D technical team, covering engineering design, project management, testing and commissioning, validation and consulting services. Together with Truking Technology, Germany ROMACO Group and Sichuan Pharmaceutical Design Institute, the company jointly undertakes EPCM general contract services for pharmaceutical enterprises, providing engineering process optimization and product lifecycle management and services for pharmaceutical enterprises.

In 2021, Truking Watertown pharmaceutical equipment Intelligent Manufacturing Industrial Park project started, located in Jiutai Economic and Technological Development Zone of Changchun Medicine Industrial Park. The new plant occupies 361,117 square meters. It is planned to complete the overall construction in three phases within 5 years. After the completion of the construction, Truking Watertown R&D Center building, intelligent workshop, employee apartment and business reception center and other modern intelligent supporting utilities will be used. The first phase of the project covers a construction area of 60,000 square meters and is expected to be put into operation in July 2023. As base of Truking Watertown with an annual output of 5,000 sets, pharmaceutical equipment intelligent manufacturing Industrial Park will greatly improve the enterprise's production capacity after the completion of the project, which will further promote the realization of the strategic objectives of the enterprise.

We insist on the core value of being respected person, making respected product and running respected enterprise and adhere to the philosophy to be the unique or to be the first and promote the spirit of Because of persistence, we are superexcellent to make Truking Watertown be one of the leading pharmaceutical water system enterprises all over the world.

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
TRUKING WATERTOWN is committed to providing compliant, efficient and reliable biological waste liquid inactivation solutions for pharmaceutical companies, which can meet the requirements of various production modes including laboratory, small-scale, pilot-scale, and large-scale production. Provide reliable support for the compliant discharge of waste liquid in the production process of the pharmaceutical industry, and promote the healthy development of the pharmaceutical industry to a green circular energy economic model.

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

The level of biological waste liquid inactivation system is divided according to the level of biosafety. Biosafety level is determined according to different degrees of biohazard, including requirements for laboratory personnel, laboratory and even environmental protection. Generally, the safety of microbial and biomedical laboratory is divided into 4 levels according to biohazard level, namely: BSL-1, BSL-2, BSL-3, BSL-4.

生物安全等级 Biological safety level	生物危害程度 Degree of biological hazard	举例 Example
BSL-1	<b>低个体危害，低群体危害 / Low individual harm, low group harm</b> 对人体、动植物或环境危害较低、不具有对健康成人、动植物致病的致病因子。 Low harm to human body, animals and plants or environment, no pathogenic factors for healthy adults, animals and plants.	<b>大肠杆菌的非致病性菌株</b> Non-pathogenic strains of E. coli
BSL-2	<b>中等个体危害，有限群体危害 / Moderate individual harm, limited group harm</b> 对人体、动植物或环境具有中等危害或具有潜在危险的致病因子、对健康成人、动植和环境不会造成严重危害。具有有效的预防和防止措施。 Disease-causing agents that are moderately harmful or potentially dangerous to human body, animals and plants or the environment, and do not cause serious harm to healthy adults, animals and the environment. With effective prevention and prevention measures.	<b>金黄色葡萄球菌</b> Staphylococcus aureus
BSL-3	<b>高个体危害，低群体危害 / High individual hazard, low group hazard</b> 对人体、动植物具有高度危害性，通过直接接触或气溶胶使人传染上严重的甚至致命疾病，或对动植物和环境具有高度危害的致病因子，通常有预防和治疗措施。 It is highly harmful to human body, animals and plants, and can infect people with serious or even fatal diseases through direct contact or aerosol, or pathogenic factors that are highly harmful to animals and plants and the environment. Prevention and treatment measures are usually available.	<b>结核分枝杆菌</b> Mycobacterium tuberculosis
BSL-4	<b>高个体危害，高群体危害 / High individual hazard, high group hazard</b> 对人体、动植物或环境具有高度危害性，通过气溶胶途径传播或者传播途径不明，或未知的、高度危险的致病因子，没有预防和防治措施。 It is highly hazardous to the human body, animals and plants or the environment. It is transmitted by aerosol route or unknown route of transmission, or unknown and highly dangerous pathogenic factors, and there is no prevention and control measures.	<b>埃博拉病毒 马尔堡病毒</b> Ebola, Marburg

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Chemical structures such as bacterial proteins and nucleic acids are connected by hydrogen bonds, and hydrogen bonds are weaker chemical bonds. When in a high temperature environment, the hydrogen bonds are destroyed, and the structures of proteins, nucleic acids, enzymes, etc. are also destroyed. Protein denaturation or coagulation, enzyme inactivation, so that bacteria and viruses lose their biological activity, resulting in death.

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For biological inactivation process, its core parameters are temperature and time, and the equipment process design is usually based on F0 numerical calculation.

F0 value is the time (min) during which the sterilization effect at a certain sterilization temperature (T) and Z at 10°C are the same as the sterilization effect at 121°C and Z at 10°C. That is, regardless of temperature changes, the effect of sterilization in t minutes is equivalent to the effect of sterilization in F0 minutes at 121 °c.

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
Running processes: biological waste liquid through a solid waste treatment unit (optional) enter the buffer unit, through the pump to heat units and heating units heated to the sterilization temperature, the heating liquid into the circulation heat preservation unit internal rules set temperature keeping time after sterilization effect, circulation cooling after discharge, waste liquid in the system continuous discharge, without stop.

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- Energy saving and consumption reduction: after the system is inactivated, the heat emitted is recovered and preheated in the system for use, resulting in high energy use efficiency.
- Uniform heating: waste liquid is heated by continuous flow, heating is uniform, there is no heating temperature gradient problem.
- Device volume: Continuously processes wastewater without storing it, eliminating the need for a large buffer tank.
- Strain capacity: can precisely adjust the temperature without changing equipment, suitable for test parameters change.
- Automatic cleaning: The system is equipped with a CIP unit that automatically descales at regular intervals.
- Flexible sterilization: block inactivation and partition inactivation are optional. Key components can be sterilized separately.
- Process parameters: equipped with a paperless recorder, you can record and verify the sterilization temperature (sterilization time) set.
- Redundant design: key components are redundant design, single pump failure can be isolated and sterilized, disassembly and maintenance after sterilization.
- Structure design: modular equipment design, the overall frame type skid-mounted structure, the bottom set overall water tray, collect accidental leakage of waste liquid equipment operation and maintenance of waste water for retreatment.

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The holding time must be greater than or equal to the processing time at the temperature corresponding to the F0 value.

For example, the treatment temperature is set to 135℃, the F0 value is 30, and the corresponding treatment time is 90s. Therefore, the heat preservation and inactivation unit needs to maintain the flow time of waste liquid ≥90s.

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
Three-dimensional pipe bending, reduce solder joints, reduce leakage risk, improve pipeline consistency.

Hard shell insulation structure, strong and durable.

Unfixed connection end, easy to maintain.

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
Single channel design reduces scaling risk. At the same time the whole removable, can be cleaned.

The two fluid media flow concentrically in spiral, the hot and cold fluid completely countercurrent heat transfer, the hot medium enters from the center hole by inward outflow, the cold fluid enters from the edge tangential flow, high heat transfer efficiency.

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The pump drive adopts magnetic drive, the pump chamber is completely isolated from the outside, and there is no sealing. Avoid the risk of waste leakage.

The redundant design of double pumps ensures continuous operation of the system.

Single pump failure can be disassembled and maintained after regional sterilization

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Continuously intercept large particle impurities from upstream waste liquid entering the inactivation system with an interception accuracy ≤1mm.


After the waste reaches the limit, steam sterilization operation is carried out, and the waste is discharged after sterilization.

Automatic operation of the system, no personnel operation for waste solid treatment, reduce personnel safety risk.

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
Biological waste liquid through a solid waste processing unit (optional) enter the buffer unit, through the pump to inactivate unit 1, after being inactivated unit 1 at level set, waste liquid into the inactivated in unit 2, inactivated heating unit 1 to sterilization temperature, maintain a set time for inactivated processing, waste liquid after entering cooling unit circulation cooling discharge after inactivated, CIP cleaning operation can be carried out in the tank after discharge. Inactivated unit 1 and inactivated unit 2 run in batches to continuously inactivate the waste liquid.

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
Principle: The injection of steam sucks the liquid from the annular hole, mixes it, and then sprays the heated hot water into the tank.  
Features: The circulation caused by the injector ensures adequate mixing and avoids temperature stratification differences.  
The sound-absorbing structure can effectively avoid the noise and vibration caused by steam injection heating.

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
- Features:
- The design of steam injection and temperature detection points for the tank volume design reduces the temperature gradient and improves the uniformity of temperature detection.
  - Optional steam tangential injection to promote waste liquid mixing.
  - The Material of Construction for the tank inner liner is S31603, and the surface is passivated.
  - Optional heating/cooling jacket design.
  - Optional mechanical mixing or magnetic mixing.
  - Spray ball is used in CIP cleaning.

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