TRUKING's Overall Solution for APIs

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Equip global pharmaceutical manufacturing industry with intelligent wings.

原料药整体解决方案 - 01 -
Overall Solution for APIs Production

全自动 CIP 系统 - 08 -
Automatic CIP system

全自动配液系统 - 09 -
Fully Automatic Preparation System

过滤洗涤干燥机 - 10 -
Filter Washing and Drying Machine

真空冷冻干燥机 - 12 -
Vacuum Freeze-dryer

冻干出料系统 - 15 -
Freeze Drying Unloading System

铝桶清洗灭菌混合分装系统 - 19 -
Aluminum Drum Cleaning, Sterilization, and Mixing-Filling System

铝桶清洗机 - 19 -
Aluminum Drum Washer

脉动真空蒸汽灭菌柜 - 21 -
Pulsating Vacuum Steam Sterilization Cabinet

单锥混合机 - 22 -
Vertical Agitator Machine

提升混合机 - 24 -
Lifting Mixer

铝桶自动转运、分装、轧盖系统 - 26 -
Aluminum drum Automatic Conveying, Filling, and Capping System

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Founded in 2000, Truking Technology Limited has become a leading pharmaceutical equipment company in China and renowned worldwide. The company specializes in pharmaceutical equipment and integrated solutions and takes the initiative to push forward the R&D of smart pharmaceutical factories. It is a listed A-share company with multiple wholly-owned or holding subsidiaries, such as Germany-based Romaco Group, Truking Watertown, Sichuan Pharmaceuticals Design Institute, Truking Feiyun, Truking Ingenuity, Truking Micro-Sphere, Truking SUT, Truking Gene, Wachine, Truking Changxing, Truking Scientific Instrument, Truking Jingbang Engineering Technology, Truking Pitide Biotechnology, Truking Boyuan Intelligent Technology, Truking New Material, Truking Intelligent Robot etc. Truking has over 9000 employees around the globe, with a total asset of over 11 billion RMB.

Truking has two operational headquarters in Changsha, China, and Germany, with research and development institutions including Changsha Central Technology Research Institute, Europe Technology Research Institute, and Sichuan Pharmaceuticals Design Institute. It has established many platforms for technology innovation, such as the National Enterprise Technology Center, National Innovation Center, Post-Doctoral Scientific Research Station, and Academician Expert Workstation. By December 31, 2023, a total of 5,466 patent applications from China have been submitted, with 4,313 patents granted and 3,090 patents currently in force. Additionally, 54 PCT international patent applications were filed, resulting in 26 patents granted in multiple countries including the United States, Russia, India, South Korea, Germany, Indonesia, Japan, and Europe. The company has led the development of over 20 national industry product technical standards for this system. The products and services of the group have now reached over 180 countries and regions, with its international market share steadily increasing year by year.

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Overall Solution for APIs Production

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Freeze dryer: it is mainly used in the pharmaceutical field such as APIs, to realize freeze-drying and drying.

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
Isolator: Autoclave unloading Isolator

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
Vacuum pump set: provide negative pressure

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CIP station: it is mainly used to clean the tank bodies and pipes of various preparation systems, freeze dryers, and reaction crystallization tanks.

无菌原料药冻干精干包整体解决方案 Overall Solution for Freeze-drying, Refining, Drying and Packaging of Aseptic APIs

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Filtering, washing and drying machine: it can function feeding, filtering, washing, drying, cooling, unloading, online cleaning and sterilization, etc. It is mainly used in the production of pharmaceutical intermediates and APIs, bio-pharmaceuticals and fine chemicals.

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
CIP cleaning station: it is mainly used to clean the tank bodies and pipes of various preparation systems, filtering, washing and drying machines, and reaction crystallization tanks.

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
Crystallization tank: it is used for material stirring and mixing, heating and cooling, finished product crystallization, etc.

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Vacuum pump set: provide negative pressure

无菌原料药溶媒结晶精干包整体解决方案 Overall Solution for Crystallization, Refining, Drying and Packaging of Aseptic APIs

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
Vacuum pump set: provide negative pressure.

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
Compounding isolator: including VHP sterilization section, vacuum feeding section, and output buffering section.

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Drum: mainly used for temporary storage, transfer and mixing of materials.

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
Aluminum drum cleaning, sterilization, transshipment and sorting system: it is mainly used in the aseptic APIs production line in pharmaceutical factories. In linked production, it can complete multiple processes such as cleaning, sterilization, material sorting, automatic sampling, and capping of aluminum drums.

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Lifting mixer: it is mainly used for drug mixing, with such functions as lifting, mixing, and station switching.

无菌原料药配料混粉分装整体解决方案 Overall Solution for Compounding, Powder Mixing and Sorting of Aseptic APIs

全自动 CIP 系统 Automatic CIP system

主要用途 / Main Application

- 主要用于各种配置系统、冻干机，过滤洗涤干燥机，各类反应结晶罐的罐体、管道清洗。
- + It is mainly used to clean the tank bodies and pipes of various configuration systems, freeze dryers, filtering, washing and drying machines, and reaction crystallization tanks.

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Fully Automatic Preparation System

应用领域 / Application

- 配液系统是按照客户的个性化要求进行设计，由理念先进的团队和专业的制造、现场施工、验证团队来实施，符合中国新版 GMP、EUGMP、FDA、cGMP 等的要求。
- 化学制药领域
 - (1) 无菌原料药配液系统
 - (2) 冻干配液系统
 - (3) CIP/SIP
- + The preparation system is designed according to the customer's individual requirements, and manufactured by an advanced concept team and a professional manufacturing, on-site construction, and verification team. It meets requirement of latest version of GMP, EUGMP, FDA and cGMP.
- + Chemical Pharmaceutical Field
 - Aseptic APIs production system
 - Freeze-drying products preparation system
 - CIP/SIP

过滤洗涤干燥机 Filter Washing and Drying Machine

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The multifunctional filter washing and drying machine completes the processes of filtration, washing, drying, and outfeed sequentially within a sealed reaction vessel. It features a simplified process flow, enhanced production efficiency, prevention of material contamination, online cleaning capability (CIP), steam sterilization (SIP), convenient material replacement, and electromechanical integration.

It finds extensive applications in sterile API, chemical, food, electronics, and other industries.

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- + Using built-in hydraulic tubing for a sleeker exterior design;
- + The filter adopts SUS316L metal woven multi-layer sintered mesh, supporting customized sintered metal mesh with a minimum pore size of 2μm; increasing filtration speed with pressure filtration or vacuum filtration;

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- Seven layers of metal sintered mesh are used, eliminating porous plates. The metal sintered mesh is connected by snap fasteners, making cleaning and replacement of the filter convenient;
- Utilizing laser automatic welding technology for welding the sintered mesh;
- The S-shaped agitator blade with rotatable lifting capability and a fully heated channel inside the blade enhances drying efficiency;
- The chassis features an automatic snap-lock structure, which automatically locks the chassis, reducing system leaks caused by differences in manual locking sequence and locking force;
- Narrower outfeed valve design and lower clearance between the agitator blade and the filter mesh result in a outfeed rate higher than industry standards;
- Optional conductivity meters and TOC probes enable online monitoring of cleaning water and organic carbon;
- A solvent crystallization turnkey solution, upstream can be matched with Truking's liquid preparation engineering for dissolution reactions, crystallization systems, and downstream can be matched with Truking's API aluminum drum cleaning, sterilization, handling, mixing, filling, and capping assembly line.

技术参数 / Technical Specification

型号 NFD Area	200	400	500	600	800	1000	1200	1400	1600	1800	2000	2300	2600
过滤面积 (m ²) Filtration Area (m ²)	0.03	0.13	0.20	0.28	0.50	0.79	1.13	1.54	2.01	2.54	3.14	4.15	5.31
总容积 (m ³) Total Volume (m ³)	0.015	0.07	0.15	0.22	0.46	0.95	1.20	1.90	2.59	3.20	4.50	6.40	9.00
搅拌电功率 (KW) Mixing Power (KW))	0.37	0.55	2.2	2.2	3	5.5	7.5	11	15	15	18.5	22	22
容器设计压力 (Mpa) Vessel Design Pressure (Mpa)	- 0.1~0.35												
装液系数 Liquid Filling Coefficient	40%~70%												
容器设计温度 (°C) Vessel Design Temperature (°C)	容器设计温度：0~150 Vessel Design Temperature: 0~150												
搅拌转速 (rpm) Mixing Speed (rpm)	0~15												
有效容积 (m ³) Effective Volume (m ³)	0.01	0.05	0.11	0.15	0.32	0.67	0.84	1.33	1.81	2.24	3.15	4.48	6.30
滤饼厚度 (mm) Filter Cake Thickness (mm)	50	120	140	180	200	220	250	300	300	350	400	400	400
搅拌桨行程 (mm) Mixing Blade Stroke (mm)	60	150	160	200	220	250	300	350	350	400	450	450	450

真空冷冻干燥机 Vacuum Freeze-dryer

主要用途 / Main Application

- ◉ 主要用于生物制品、原料药（天然药物、热敏性药物、抗生素等）等制药领域，实现冻干过程，达到干燥目的。
- + It is mainly used in pharmaceutical fields such as biological products, APIs (natural medicines, heat-sensitive medicines, antibiotics, etc.), to realize the freeze-drying and drying.

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
Membrane freeze-drying

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Stainless steel tray loading and unloading


技术参数 / Technical specification

面积 Area	1	3	5	10	15	20	30	40
有效面积 (m ²) Effective area (m ²)	1.00	2.88	5.00	9.51	15.80	20.21	30.10	40.12
隔板数量 (个) Number of separators (pcs)	3+1	4+1	6+1	6+1	10+1	13+1	12+1	16+1
隔板温度范围 (°C) Separator temperature range (°C)	(-55°C~ +8°C)							
冷凝器温度 (°C) Condenser temperature (°C)	≤ -75°C							
极限真空 (mbar) Ultimate vacuum (mbar)	0.005							
系统泄漏率 (mbar.L/s) System leakage rate (mbar.L/s)	0.01							

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Notes: Special size can be customized.

托盘冻干出料系统 Tray Freeze-Drying Unloading System

主要用途 / Main Application

- 托盘原料药冻干出料系统，适用于高端生物制剂的生产。
- + Tray APIs Freeze-Drying Unloading System, suitable for the production of high-end biopharmaceuticals.

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
Flanged-Shelf Freeze-Drying Unloading System

主要用途 / Main Application

- ◉ Orabs 原料药全自动真空吸料输送出料系统，适用于规模化大生产。自动出料可一套系统对应对多台冻干机。
- ◉ Isolator 原料药全自动真空吸料输送出料系统，适用于致敏性、带毒及含水率要求较高的产品
- + Orabs APIs Fully Automatic Vacuum Suction Feeding and Unloading System, suitable for large-scale production. One automatic unloading system can correspond to multiple freeze dryers.
- + Isolator API Fully Automatic Vacuum Suction Feeding and Unloading System, suitable for sensitizing, toxic, and high-moisture-requirement products.

性能特点 / Features

- ◉ 大规模生产，全自动解决方案；
- ◉ 可实现一对多处理；
- ◉ 开小门等高出料，配合充氮处理，可有效控制成品含水率。
- + Mass production with a fully automated solution;
- + Enables one-to-many processing;
- + Features uniform height unloading with pizza door openings, combined with nitrogen processing to effectively control the moisture content of the final product.

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The Aluminum Drum Washer is primarily used for cleaning raw aluminum drums.

The Truking Aluminum Drum Washer is a fully automatic cleaning device, consisting of automatic barrel feeding rollers, an ultrasonic cleaning tank, purified water rinse + compressed air blowing, injection water cleaning, and compressed air blowing. The cleaning process requires no manual intervention, thus reducing the risk of human error.

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脉动真空蒸汽灭菌柜 Pulsating Vacuum Steam Sterilization Cabinet

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
Pulsating vacuum steam sterilization cabinet is a machine integrating vacuum technology and fluid technology. It is composed of box, door, frame, outer sheet metal, vacuum system, steam inlet system, exhaust system, compressed air system, transfer system, filter online sterilization system, hot air system, control system and other functional branch systems. The pulsating vacuum steam sterilization process can effectively remove the non condensable gas in the chamber, so as to ensure that the steam can penetrate into the porous loose load quickly and effectively, and ensure that the temperature in the cabinet is uniform. After sterilization, the product is dried by vacuumizing and adding clean and dry air to remove water vapor in the cavity. Process flow: preheating - pulsating - heating - sterilization - drying and cooling

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The Truking VAM type single cone mixer integrates mixing and drying functions into one unit. It enables vacuum feeding, online CIP (Cleaning In Place), SIP (Sterilization In Place), and can be paired with gravity or vibration dosing systems for filling.

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The main purpose of this equipment is for drug mixing in pharmaceutical factories, with functions including lifting, mixing, and workstation switching.

A lifting mixer can be equipped with mixing hoppers of different specifications to meet the mixing requirements of large quantities and multiple varieties.

性能特点 / Features

- 电动升降、料桶电动混合、底盘电动旋转，采用全伺服系统，相比液压站更干净整洁；
 - 易清洁，无泄漏，适合于无菌室；
 - 升降采用电动丝杆驱动，无液压油的泄露而带来的料桶沉降的风险；
 - 升降配置丝杆和抱闸伺服电机；
 - 混合电机采用伺服电机，可设置正反转时间及转速，达到最佳的混合效果。
- + Electric lifting, electric mixing of barrels, and electric rotation of the chassis are achieved through a fully servo system, which is cleaner and tidier compared to hydraulic stations;
 - + Easy to clean, leak-free, suitable for sterile environments;
 - + The lifting mechanism is driven by an electric screw rod, eliminating the risk of barrel settling caused by hydraulic oil leakage;
 - + The lifting system is equipped with a screw rod and brake servo motor;
 - + The mixing motor employs a servo motor, allowing for adjustable mixing time and speed, achieving optimal mixing effects.

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Aluminum Drum Automatic Conveying, Filling, and Capping System

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Primarily designed for use in API product lines, this system seamlessly interfaces with various processes including aluminum drum cleaning machines, sterilization cabinet robotic arm, filling systems, capping, and more.

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冻干工艺服务 Freeze-drying Process Service

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冻干显微镜
Freeze-drying Microscope

差示扫描量热仪
Differential scanning calorimeter

- Truking freeze-drying laboratory conducts in-depth research on the entire process of freeze-drying of products. It provides the best freeze-drying process solution which meets the requirements of GMP, and interprets the relationship between the production quality, efficiency and regulations.

实验型冻干机
Laboratory Freeze Dryer

提供的服务 / Service can be Provided:

- 工艺开发：帮助客户开发冻干曲线；
 - 处方设计：帮助客户开发药品冻干的处方；
 - 工艺放大：为需要从小规模生产或实验小试转变为大规模生产的客户提供帮助；
 - 工艺优化：帮助已拥有原始工艺曲线且需要进一步完善的客户，优化工艺曲线，提升药品冻干速率，缩短冻干周期。
- Process development: Help user to develop freeze-drying curve;
 - Recipe design: Help user to develop recipe of freeze-drying product;
 - Process scale up: Offer help to user who need to transform from small scale production or laboratory to large scale production;
 - Process optimization: Help user who has original process curve, but need to optimize, to optimize the process curve, improve freeze-drying speed, and reduce freeze-drying cycle.

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Service Concept

400-9988-900

Adhering to concept of “all for customers, for all customers” , and with the mission of “creating value for customers” , TRUKING focuses on building first-class R&D, design, sales, manufacturing, service teams and service management system to offer high-efficient and high-quality products and services for the customers.

◎ Project Planning Support

TRUKING has a series of professional medical and pharmaceutical design institutes, product research institutes, a group of senior scientific research and technical personnel, national-level aseptic pharmaceutical process verification experts, professors, etc., to provide scientific and effective technical planning support for customers' projects.

◎ Process Assurance Service

TRUKING launched the overall solution for medical & pharmaceutical equipment realizing the full automation and aseptic production process, intelligent detection and sorting of the drug, as well as computer verification assurance so as to produce assured products for market.

◎ Service In The Park

As an equipment supplier, TRUKING not only integrates industrial regulations and standards into product design and manufacturing, but also introduces the industrial regulations and standards into Truking industrial park in the form of lectures and exchanges to interact and integrate with the technology and production of pharmaceutical enterprises.

◎ Equipment Management Consultation Service

TRUKING aims to satisfy the maximum benefited demands of the customer and to provide a package of equipment management services to realize minimization investment and maximization production output together with pharmaceutical factories.

◎ Adhere to the philosophy of "customer-centric", TRUKING creates the extraordinary service experience with the spirit of the craftsman, to create more value for customers.

◎ We are always ready to escort the normal operation of the customer's equipment.

- + Onsite equipment installation, commissioning, acceptance
- + Training for equipment process, operation and maintenance
- + Periodic door-to-door maintenance and return visit
- + Performance test and maintenance for machinery parts
- + Timely supply of spare parts
- + Product technology upgrading and secondary development
- + Equipment overhaul

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