**Truking Intelligent
sterilization/cleaning technology**

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Equip global pharmaceutical manufacturing industry with intelligent wings.

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Founded in 2000, Truking Technology Limited has become a leading pharmaceutical equipment company in China and renowned worldwide. The company specializes in pharmaceutical equipment and integrated solutions and takes the initiative to push forward the R&D of smart pharmaceutical factories. It is a listed A-share company with multiple wholly-owned or holding subsidiaries, such as Germany-based Romaco Group, Truking Watertown, Sichuan Pharmaceuticals Design Institute, Truking Feiyun, Truking Ingenuity, Truking Micro-Sphere, Truking SUT, Truking Gene, Wachine, Truking Changxing, Truking Scientific Instrument, Truking Jingbang Engineering Technology, Truking Pitide Biotechnology, Truking Boyuan Intelligent Technology, Truking New Material, Truking Intelligent Robot etc. Truking has over 9000 employees around the globe, with a total asset of over 11 billion RMB.

Truking has two operational headquarters in Changsha, China, and Germany, with research and development institutions including Changsha Central Technology Research Institute, Europe Technology Research Institute, and Sichuan Pharmaceuticals Design Institute. It has established many platforms for technology innovation, such as the National Enterprise Technology Center, National Innovation Center, Post-Doctoral Scientific Research Station, and Academician Expert Workstation. By December 31, 2023, a total of 5,466 patent applications from China have been submitted, with 4,313 patents granted and 3,090 patents currently in force. Additionally, 54 PCT international patent applications were filed, resulting in 26 patents granted in multiple countries including the United States, Russia, India, South Korea, Germany, Indonesia, Japan, and Europe. The company has led the development of over 20 national industry product technical standards for this system. The products and services of the group have now reached over 180 countries and regions, with its international market share steadily increasing year by year.

脉动真空蒸汽灭菌柜 Pulsating Vacuum Steam Sterilizer

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This series of equipment adopts saturated pure steam as the sterilization medium, to do sterilization and drying treatment on clothing, dressings, instruments, culture medium, glass, etc. The equipment strictly follows the GMP requirements during the process of design, manufacture, test and acceptance, and also online detection and real-time printing of temperature and pressure in the whole process of sterilization, the whole process can be traced. It is widely used in hospitals, pharmaceutical factory, laboratory, scientific research units and so on for disinfectant and sterilization.

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+ Unique sealing design:

Door sealing: structure of dovetail groove bi-directional sealing design, to ensure seamless connection during switch between pulse vacuum and steam sterilization, avoiding leakage risk of sterilization chamber.

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+ Seal-off Air-Pressure Maintenance System: Equipped with sealing buffers and automatic control systems, it achieves a secure self-sealing pressure maintenance even after equipment shutdown and air pressure cutoff, ensuring the safety of personnel, equipment, and products.

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+ Pure vapor diffusion panel design of the chamber, can control condensed water distributary and steam sterilizing quality. And coupled with extensive testing and validation of hydrophobic throttle orifice plate joint design, with higher sterilization efficiency and better effect.

+ TIG welding is applied for the chamber and pipeline, mirror polishing to $Ra \leq 0.5\mu m$.

+ Flange inside the chamber adopts overall processing technology, without welding residual stress. Advanced control of the sealing surface and sealing groove size during the manufacturing process, assures greater safety and reliability of the sealing during operation, and thus service life of the sealing ring is lengthened.

+ Automatic detection control of pressure and temperature inside chamber, Sealing door and the whole inner cavity with insulation structure, to improve temperature distribution.

+ High density ceramic insulation material is used, no particle dropping. And covering Stainless steel wire drawing board outside.

+ It can realize SIP and integrity test for sterilizing filter, providing greater safety and reliability during usage.

+ Air Detection System -ADS: Based on the sterilization program in operation, it automatically captures the mixed gas inside the sterilization equipment, determines through an intelligent system whether the content of non-condensable gas in the captured gas affects the sterilization effect, and generates electronic batch records.

+ Real-time monitoring of temperature and pressure, and the curve is real-time printed and recorded.

+ Portal automatic control system completely meets the requirements of GAMP5 and 21CFRPart11.

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This type of equipment adopts mixture of air and steam as heating media, and is particularly suitable for sterilization of products with easy-deformed packing. It strictly follows the GMP requirements for design, manufacture, testing and acceptance, which can realize on-line monitoring and real-time print of temperature and pressure in sterilization process, the whole process can be traced. Typically used for terminal sterilization of BFS and prefilled syringes products.

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Unique sealing design:

Door sealing: structure of dovetail groove bi-directional sealing design to ensure seamless docking at all stages of sterilization, avoiding leakage risk of sterilization chamber.

Sealing system: pressure protection of inflatable sealing system, with pressure detection alarm function.

Outlet steam separator: adopt air partition structure and with steam pressure relief outlet to prevent pollution risks brought by recoil of the floor drain pipeline.

Main circulation fan system is driven by magnetic structure, can ensure the sterilization chamber is absolute separated and completely sealed, with easy maintenance, low working and maintenance costs, and smooth operation, and it has low noise, stable performance and long service life.

Chamber design has separate water cooling system, with fast cooling. Cooling capacity can be adjusted through the PID Program, and with a pressure protection function. Cooling water can be recycled for saving energy and water. Pipe and valves in cooling system are easily installed and dismantled, is convenient for regular cleaning for the inner chamber of equipment and pipe cleaning for the cooling water pipeline.

Pure vapor diffusion panel design of the chamber, can control condensed water tributary and steam sterilizing quality. And coupled with extensive testing and validation of hydrophobic throttle orifice plate joint design, with higher sterilization efficiency the better effect.

TIG welding is applied for the chamber and pipeline, mirror polishing to $Ra \leq 0.5\mu m$.

Chamber flange adopts overall processing technology, without Welding residual stress. Advanced control of the sealing surface and sealing groove size during the manufacturing process, assures greater safety and reliability of the sealing during operation, and thus service life of the sealing ring is lengthened.

Automatic detection and control of the pressure and temperature inside the chamber, Sealing door and the whole inner cavity with insulation structure, to improve temperature distribution. High density ceramic insulation material is used, no particle dropping. And Stainless steel wire drawing board outside.

It can realize SIP and integrity test for sterilizing filter, providing greater safety and reliability during usage.

Real-time monitoring of temperature and pressure, and the curve is real-time printed and recorded.

Portal automatic control system completely meets the requirements of GAMP5 and 21CFRPart11.

产品参数 / Product Parameters

主要技术参数 Main Parameters			
设计压力 Design Pressure	0.3/-0.1MPa	纯蒸汽源压力 Pure steam source pressure	0.2-0.4MPa
工作压力 Working Pressure	0.25/-0.097MPa	工业蒸汽源压力 Industrial steam source pressure	0.4-0.6MPa
设计温度 Design Temperature	144°C	水源压力 Water pressure	0.2-0.4MPa
工作温度 Working Temperature	139/20°C	压缩空气压力 Compressed Air Pressure	0.6-0.8MPa

规格型号 / Specification model

参数 Parameters	容积 Volume	净重 Net weight	内室尺寸 Room size	外形尺寸 Dimensions	纯蒸汽耗量 Pure steam consumption	工业蒸汽耗量 Industrial steam consumption	水耗量 Water consumption	压缩空气耗量 Compressed air consumption	总电源 Power
型号 Model	m³	Kg	WxHxD (mm)	WxHxD (mm)	Kg/cycle	Kg/cycle	Kg/cycle	m³/h	Kw
SAS0405	0.4	1950	520x760x1000	1735x1970x1300	15	20	500	12	6
SAS0505	0.5	2200	520x760x1250	1735x1970x1550	16	22	560	12	6
SAS0706	0.7	2400	640x1010x1050	1855x2350x1350	18	25	640	12	6
SAS0806	0.8	2600	640x1010x1250	1855x2350x1550	20	28	700	15	6
SAS1008	1	2800	840x1010x1250	2055x2350x1550	25	30	780	15	6
SAS1408	1.4	2950	840x1010x1700	2055x2350x2000	28	35	860	25	6
SAS1808	1.8	3100	840x1010x2100	2055x2350x2400	32	38	1150	25	10
SAS2108	2.1	3500	840x1010x2500	2055x2350x2800	40	45	1280	28	10
SAS2608	2.6	4000	840x1010x3100	2055x2350x3400	48	52	1850	28	12
SAS3108	3.1	4500	840x1010x3700	2055x2350x4000	60	65	2200	35	13

水浴灭菌柜 Super-heated Water Sterilizer

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This series of equipment adopts super-heated circulating water as the sterilization medium, to do spray or shower sterilization on filled product. (Such as ampoules, oral liquid bottles, glass bottles, plastic bottles, soft bags etc.) The equipment strictly follows the GMP requirements during the process of design, manufacture, test and acceptance, and also realizes online monitoring and real-time printing of temperature and pressure in the whole process of sterilization, the whole process can be traced. It is widely used in pharmaceutical factory, laboratory, scientific research units and so on to realize disinfectant and sterilization.

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+ Unique sealing design:

Door sealing: structure of dovetail groove bi-directional sealing design to ensure seamless connection during switch between pulse vacuum and steam sterilization, avoiding leakage risk of sterilization chamber.

Sealing system: pressure protection of inflatable sealing system, with pressure detection alarm function

+ Outlet steam separator: adopt air partition structure and with steam pressure relief outlet to prevent pollution risks brought by recoil of the floor drain pipeline.

+ The main body is analyzed by finite element, achieving the optimal structural strength. The compressed air diffusion panel design of the chamber, can control distributaries of compressed air, reduce the chamber temperature and pressure fluctuation, with higher sterilization efficiency and better effect.

+ TIG welding is applied for the chamber and pipeline, mirror polishing to $Ra \leq 0.6\mu m$

+ Chamber flange adopts overall processing technology, without welding residual stress. Advanced control of the sealing surface and sealing groove size during the manufacturing process, assures greater safety and reliability of the sealing during operation, and thus service life of the sealing ring is lengthened.

+ Automatic detection and control of the pressure and temperature inside the chamber, Sealing door and the whole inner cavity with insulation structure, to improve the temperature distribution. High density ceramic insulation material is used, no particle dropping. And Stainless steel wire drawing board outside.

+ Hot water spray sterilization process is the most effective and economic process.

+ Sterilizer chamber, double door, heat circulation pipe are designed and manufactured according to the standard of pressure vessel.

+ It has two systems: manual transfer and automatic transfer, with stable running, high efficiency and less operators.

+ The control cabinet has display system, provides real-time operation status, with special port for verification.

+ When the equipment is used for sterilization of lipid emulsion infusion or suspension infusion, it is equipped with interior rotary device in the chamber, with reliable sealing for shaft hole.

+ Controlled by industrial computer, with professionally designed program. It can provide different sterilization program as per user's request, it has F0 value processing system, and temperature

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pressure, F0 value and their curve graph are constantly monitored and recorded in sterilization, and displayed in real-time.

- + Tracking control in the process of heating and cooling, higher than the preset and the fixed pressure difference of saturated vapor pressure, maintain via compressed air, to prevent breakage of glass bottle and deformation of soft packing.
- + For the sterilization of emulsion products, a rotary sterilizer is configured to control the parking of the rolling body in situ to ensure that the transportation track of the sterilization vehicle inside and outside the cabinet is smooth and the sterilization vehicle enters and exits the cabinet smoothly;
- + Real-time monitoring of temperature and pressure, and real-time printing and recording of curve.
- + Portal automatic control system completely meets the requirements of GAMP5 and 21CFRPart11.

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- + Safety devices are installed on the door, automatically stop when encountering obstacles. Door is equipped with sealed pressure balance device, can remain sealing pressure in case of circuit fault;
- + Door is equipped with safety interlock, to prevent the door open in the following circumstances: (1) pressure in chamber; (2) sterilization temperature higher than 60 °C ; (3) level of circulation water is higher than the water collecting sump.
- + Inlets of steam, circulation water, cooling water, compressed air pipeline are equipped with filters, ensure safety condition for executing components.

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PW series Parts cleaner is a stable and efficient cleaning equipment designed, manufactured, tested and accepted in strict compliance with cGMP, GAMP, FDA and EU standards; it is suitable for: glassware, plastic and stainless steel containers, production line tools such as Cleaning and drying of various materials such as pipes and rubber plug hoppers. Widely used in bioengineering, animal testing, pharmaceutical industry and other fields;

The automatic cleaning machine for cleaning process parameters complete monitoring, with detailed batch records, can achieve the entire cleaning step data recording, traceability and repeatability.

优势特点 / Performance advantages

- ◎ 清洗舱内表面镜面抛光至 $Ra \leq 0.4\mu m$ ，水槽及清洗舱顶采用十字筋结构设计，清洗舱底部残留水可自排尽，顶部水减少附着量，确保设备良好的自清洁效果；
- ◎ 清洗舱采用无颗粒脱落环保的保温隔热层，腔室外表面不超过环境温度 $20^{\circ}C$ ，节能减耗，保温层外覆优质 S30408 不锈钢拉丝保温板，美观且便于清洁；
- ◎ 根据需求匹配相应的卫生级循环水泵，循环泵出口带压力监测；变频控制循环水泵，根据不同的清洗需求可设定不同的喷淋水压，实现器具的高效清洁；
- ◎ 热风排放设计有除湿冷却系统，可实现汽水分离及热风排放冷却功能，降低排至夹层风管的气体中水分含量，降低夹层风管内产生积水，同时避免热风对风管内元器件腐蚀；
- ◎ 多路清洗剂注入系统，清洗剂输送管路单向截断功能，配置流量监测，通过时间监测注入液量；并配置高电导率监测，保证清洗剂注入量精确稳定。清洗剂注入数据的可追溯以及重复；
- ◎ 低电导率日常监测清洗终洗水电导率值，并进行终洗步段进行控制，确保清洗效果满足要求；
- ◎ 独立的双通道热风干燥系统，均配置高压鼓风机、加热系统、HEPA 过滤器，一路将高温洁净热风送入清洗舱内部干燥，一路送入管道并通过对接接头从清洗架喷淋球喷出，实现管路、清洗舱、清洗架及清洗物品快速干燥；
- ◎ 密封门为双层中空钢化玻璃门板，周边采用不锈钢板包边保护，且具备隔热降噪功能；
- ◎ 定制化清洗架，一对一清洗工位设计，实现对需清洁器具 100% 清洗覆盖，插销式喷嘴便于拆卸维护及清洁；
- ◎ 采用工业计算机精密控制，专业程序式设计。根据不同的清洗器具，建立不同的清洗机制，并可自由搭配不同的清洗工艺；
- + The mirror surface of the cleaning cabin is polished to $Ra \leq 0.4\mu m$. The tank and the roof of the washing cabin are designed with a cross-rib structure. The residual water at the bottom of the washing cabin can be drained away automatically, and the amount of water on the top of the cabin can be reduced to ensure a good self-cleaning effect of the equipment
- + The cleaning cabin adopts an environmentally friendly thermal insulation layer with no particles falling off. The outer surface of the chamber does not exceed the ambient temperature of $20^{\circ}C$. Energy saving and consumption reduction. The thermal insulation layer is covered with high-quality S30408 stainless steel wire drawing insulation board, which is beautiful and easy to clean;
- + Match the corresponding sanitary-grade circulating water pump according to the demand, with pressure monitoring at the outlet of the circulating pump; variable frequency control circulating water pump, according to different cleaning requirements, different spray water pressure can be set to achieve efficient cleaning of the appliance;
- + The hot air discharge is designed with a dehumidification cooling system, which can realize the functions of steam-water separation and hot air discharge cooling, reduce the moisture content of the gas discharged into the interlayer air duct, reduce the water accumulation in the interlayer air duct, and avoid the hot air corroding the components in the air duct
- + Multi-channel cleaning agent injection system, one-way shut-off function of cleaning agent delivery line, configuration of flow monitoring, monitoring of injection volume through time; and configuration of high conductivity monitoring to ensure accurate and stable injection of cleaning agent. Traceability and duplication of cleaning agent injection data;
- + Daily monitoring of the conductivity value of the final washing water with low conductivity, and control the final washing step to ensure that the cleaning effect meets the requirements.
- + The independent two-channel hot air drying system, all equipped with high pressure blower, heating system, HEPA filter, all the way to the high temperature clean hot air into the cleaning tank internal drying, all the way into the pipe and through the butt joint from the cleaning rack spray ball, to achieve the pipeline, cleaning tank, cleaning rack and cleaning items quickly drying;
- + The sealed door is a double-layer hollow tempered glass door panel, the periphery is protected by a stainless steel plate, and it has the function of heat insulation and noise reduction;
- + Customized cleaning rack, one-to-one cleaning station design, to achieve 100% cleaning coverage of cleaning parts, plug-type nozzles are convenient for disassembly maintenance and cleaning;

- 强大的数据存储能力，工业计算机至少可存储 1000 批生产批次数据，运行微软 SQL 数据库管理软件，数据安全、可靠，存储效率高；
 - 实时生成美观、数据信息丰富完善的电子数据报表和曲线报表，可以实现数据的在线打印和离线报表打印；
 - 自控系统完全符合 GAMP5 和 21CFRPart11 要求。
- + Adopt industrial computer precision control, professional program design. According to different cleaning appliances, establish different cleaning mechanisms, and can freely match different cleaning
 - + Powerful data storage capabilities, industrial computers can store at least 1,000 batches of production data, run Microsoft SQL database management software, data security, reliability, and high storage efficiency;
 - + Real-time generation of beautiful electronic data reports and curve reports with rich and complete data information can realize online printing of data and offline report printing;
 - + Portal automatic control system completely meets the requirements of GAMP5 and 21CFRPart11.

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It is mainly used for the washing of outer wall residue of vial, ampoule and oral liquid bottles. It can effectively wash outer wall residues of high-activity/highly sensitizing material (freeze-drying) powder product and products with high requirement for appearance quality to improve the performance and efficiency of light inspection and labeling.

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- + It adopts active protection of aluminum cap, preventing washing water from contaminating the gap of aluminum cap and vial body.
- + Ring nozzle achieves 360° continuously tracking washing and drying; the washing ring goes up and down around the vial body to rinse

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and blow; the bottom is separately washed and blown. The washing and blow-drying area receive good even effects.

- + It remains simple to change the specification, and different specifications enable adjustment of overall height for washing. The device adopts big turntable suspending in the air clamp structure to wash and dry for comprehensive washing and drying.
- + It employs compressed air to blow the vial body and bottom to be dry, there is no hot air drying; this ensures that the temperature of the vial does not rise so that the drug will not be affected as well.
- + It is optional that brush washing function can be configured, the bottle rotates due to the belt, silicone soft hair brushes and washes the vial body and bottom, thoroughly washing the surface of the vial body.
- + Multi bottom rail is adopted in the bottom of the product to reduce the contact area of the product and bottom rail.
- + The pipe is designed with functions of emptying compressed air and water discharge at low point.
- + The device adopts servo plus timing belt structure to reduce energy consumption with low conveyance noise.

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Adhering to concept of “all for customers, for all customers” , and with the mission of “creating value for customers” , TRUKING focuses on building first-class R&D, design, sales, manufacturing, service teams and service management system to offer high-efficient and high-quality products and services for the customers.

◎ Project Planning Support

TRUKING has a series of professional medical and pharmaceutical design institutes, product research institutes, a group of senior scientific research and technical personnel, aseptic pharmaceutical process verification experts, professors, etc., to provide scientific and effective technical planning support for customers' projects.

◎ Process Assurance Service

TRUKING launched the overall solution for medical & pharmaceutical equipment realizing the full automation and aseptic production process, intelligent detection and sorting of the drug, as well as computer verification assurance so as to produce assured products for market.

◎ Service In The Park

As an equipment supplier, TRUKING not only integrates industrial regulations and standards into product design and manufacturing, but also introduces the industrial regulations and standards into Truking industrial park in the form of lectures and exchanges to interact and integrate with the technology and production of pharmaceutical enterprises.

◎ Equipment Management Consultation Service

TRUKING aims to satisfy the maximum benefited demands of the customer and to provide a package of equipment management services to realize minimization investment and maximization production output together with pharmaceutical factories.

◎ Adhere to the philosophy of "customer-centric", TRUKING creates the extraordinary service experience with the spirit of the craftsman, to create more value for customers.

◎ We are always ready to escort the normal operation of the customer's equipment.

- + Onsite equipment installation, commissioning, acceptance
- + Training for equipment process, operation and maintenance
- + Periodic door-to-door maintenance and return visit
- + Performance test and maintenance for machinery parts
- + Timely supply of spare parts
- + Product technology upgrading and secondary development
- + Equipment overhaul

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