Truking Rhythmic R5 Smart Stick Packing Line

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Equip global pharmaceutical manufacturing industry with intelligent wings.

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Truking's Rhythmic R5 Smart Stick Packing Line

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SP-G12, 12-column Stick Packing Machine for Granules

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SP-P12, 12-column Stick Packing Machine for Powders

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SP-L12, 12-column Stick Packing Machine for Powders

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SP-W12, Multi-column stick weighing machine

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
Founded in 2000, Truking Technology Limited has become a leading pharmaceutical equipment company in China and renowned worldwide. The company specializes in pharmaceutical equipment and integrated solutions and takes the initiative to push forward the R&D of smart pharmaceutical factories. It is a listed A-share company with multiple wholly-owned or holding subsidiaries, such as Germany-based Romaco Group, Truking Watertown, Sichuan Pharmaceuticals Design Institute, Truking Feiyun, Truking Ingenuity, Truking Micro-Sphere, Truking SUT, Truking Gene, Wachine, Truking Changxing, Truking Scientific Instrument, Truking Jingbang Engineering Technology, Truking Pitide Biotechnology, Truking Boyuan Intelligent Technology, Truking New Material, Truking Intelligent Robot etc. Truking has over 9000 employees around the globe, with a total asset of over 11 billion RMB.

Truking has two operational headquarters in Changsha, China, and Germany, with research and development institutions including Changsha Central Technology Research Institute, Europe Technology Research Institute, and Sichuan Pharmaceuticals Design Institute. It has established many platforms for technology innovation, such as the National Enterprise Technology Center, National Innovation Center, Post-Doctoral Scientific Research Station, and Academician Expert Workstation. By December 31, 2023, a total of 5,466 patent applications from China have been submitted, with 4,313 patents granted and 3,090 patents currently in force. Additionally, 54 PCT international patent applications were filed, resulting in 26 patents granted in multiple countries including the United States, Russia, India, South Korea, Germany, Indonesia, Japan, and Europe. The company has led the development of over 20 national industry product technical standards for this system. The products and services of the group have now reached over 180 countries and regions, with its international market share steadily increasing year by year.

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This production line is applicable to quantitative stick packing of various solid substances, such as Powder, liquid, Granules, pill, ointmentXXXXXX. It can be connected with the stick package weighing, sorting and secondary counting devices. Application of cartoner and downstream packing line is allowed to facilitate box/container packing, palletizing of stick package on demand. It is widely used in the pharmaceutical, healthcare, food, daily chemical, milk powder and other sectors.

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条包机优势特点 / Advantages & Features

- ◎ 稳定高效
楚天强大的制造能力，90% 零件自己加工，保证加工精度。
 - ◎ 自动纠偏
整体伺服纠偏，纠偏控制更精准。膜分切精度偏差控制在 0.5mm 以内。
 - ◎ 全伺服控制
开卷、膜纠偏、色标调节、灌装、纵封、横封、拉膜牵引、圆角切，均采用伺服控制。袋长平均
- + Stable And Efficient
Truking boasts excellent manufacturing capability, who can independently complete machining of 90% parts to specified accuracy
 - + Auto Skew Adjustment
Adopt overall servo to correct deviation, deviation correction is more accurate. The film slitting accuracy deviation is controlled within 0.5mm.
 - + Full servo control
Servo control applies to uncoiling, membrane misalignment correction,color code, filling, longitudinal sealing, transverse

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sealing, avoidance at specified locations, horizontal sealing membrane-drawing and fillet cutting. The average deviation of the bag length is controlled within 0.5mm.

+ linear guide rail

The slide guide moves across the straight line guide rail, featuring more advantages in transmission accuracy, loading, stability and durability than the linear bearing.

+ high precision servo filling

Excellent machine processing level of Truking, each column adopts independent servo drive, real-time communication with multi-column weighing, filling volume of each column can be separately adjusted.

+ Easy to open at one attempt, good user experience, low material cost

Horizontal and longitudinal easy openings are formed in one mold to ensure that double easy openings stay at one straight line, thus improving user experience and ensuring they are easy to open. The average deviation of the bag length is controlled within 0.5mm.

+ Easy to replace and clean

Adopt toolless quick disassembly design. The disassembly process shall not exceed 15 minutes. The mechanism of Longitudinal sealing and horizontal sealing can be opened, easy to clean.


+ The film can be changed quickly and accurately

The film can be changed quickly and accurately Vacuum negative pressure film platform. Controlled film tiling, Equipped with guide cutter to ensure the straightness of cutting, to ensure film is flat.

+ Easy to operate and maintain.

Equipped with sky blue protective door, can be observed in all directions, no dead Angle; Large internal space, easy operation, easy maintenance.

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This production line is applicable to quantitative stick packaging of various granular substances. It is widely used in the pharmaceutical, healthcare, food and daily chemical areas etc.

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外形尺寸 Overall dimension	2800mm(L)×2550mm(W)×4300mm(H)
设备重量 Weight	2600kg
适应规格 Applicable specification	长度≤ 200mm, 宽度≤ 60mm, 袋长在 50-200mm 范围内可调 Length ≤ 200mm, width ≤ 60mm, adjustable length within 50-200mm
最高产能 Max. capacity	70 切 / 分钟 70 cutting/min
条包列数 Column no.s of stick packing	12 列 12-column
最大总膜宽 Max. total membrane width	1000mm
横封 Horizontal sealing	伺服量杯灌装 Servo volume filling
设备噪音 Noise	≤ 75 分贝 ≤ 75dB
主电源 Main power supply	三相 AC380V-50Hz-5 线 3-phase AC380V -50Hz-5-wire
额定功率 Rated power	35.0kW
压缩空气耗量 Consumption of compressed air	快插接口：Φ16mm 气管 Quick-change adaptor: Φ16mm gas hose 300L/min (0.6MPa-0.8MPa)

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This production line is applicable to quantitative stick packaging of various powder substances. It is widely used in the pharmaceutical, healthcare, food, daily chemical areas and other sectors.

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外形尺寸 Overall dimension	2800mm(L)×2550mm(W)×4300mm(H)
设备重量 Weight	2600kg
适应规格 Applicable specification	长度≤ 200mm, 宽度≤ 60mm, 袋长在 50-200mm 范围内可调 Length ≤ 200mm, width ≤ 60mm, adjustable length within 50-200mm
最高产能 Max. capacity	70 切 / 分钟 70 cutting/min.
条包列数 Column no.s of stick packing	12 列 12-column
最大总膜宽 Max. total membrane width	1000mm
横封 Horizontal sealing	伺服螺杆灌装 Servo screw filling
设备噪音 Noise	≤ 75 分贝 ≤ 75dB
主电源 Main power supply	三相 AC380V-50Hz-5 线 3-phase AC380V -50Hz-5-wire
额定功率 Rated power	35.0kW
压缩空气耗量 Consumption of compressed air	快插接口: Φ16mm 气管 Quick-change adaptor: φ16mm gas hose 300L/min (0.6MPa-0.8MPa)

SP-L12, 12 列液体条包机 SP-L12, 12-column Stick Packing Machine for Powders

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Suitable for a variety of oral liquid, mouthwash, paste zi and other substances of the quantitative strip bag packaging, can be widely used in pharmaceutical, health care, food, daily chemical and other industries.

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外形尺寸 Overall dimension	2500mm(L)×2350mm(W)×3650mm(H)- 主机 2500mm(L)×1200mm(W)×2500mm(H)- 灌装
设备重量 Weight	3000kg
适应规格 Applicable specification	长度≤ 160mm, 宽度≤ 60mm, 袋长在 80-160mm 范围内可调 Length ≤ 160mm, width ≤ 60mm, adjustable length within 80-160mm
最高产能 Max. capacity	55 切 / 分钟 55cut/min
条包列数 Column no.s of stick packing	12 列 12-column
最大总膜宽 Max. total membrane width	1000mm
横封 Horizontal sealing	双横封 Double transverse sealing
易撕口形式 The form of easy to tear	横封纵封易撕口一刀切成 Easy tear of horizontal seal and vertical seal can be cut in one time
切刀形式 The form of cutter	圆角切刀 Radius tool
灌装形式 Filling type	伺服陶瓷泵灌装 filling by Servo-ceramic pump
主电源 Main power supply	3 相 ~+ 零线 + 接地 AC380V/50Hz/60Hz 3 Phase ~ + zero line + ground AC380V / 50Hz / 60Hz
额定功率 Rated power	35.0kW
压缩空气耗量 Consumption of compressed air	快插接口: Φ16mm 气管 Quick-change adaptor: φ16mm gas hose 800L/min (0.6MPa-0.8MPa)
氮气 (有充氮要求时) Nitrogen gas (when nitrogen is required)	纯度 >99.99%, 压力: 0.4Mpa; 60 立方米 / 小时 / 台 Purity > 99.99%, pressure: 0.4Mpa; 60 m ³ / h / unit

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This equipment is used with multiple lines of charter flight wiring. Complete the orderly delivery of the bags, enter the weighing system for single column weighing detection and elimination.

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Ensure that the bag will not overrush, and the transmission is reliable and orderly.
 - + sample

The automatic sampling mechanism can choose regular sampling or manual random sampling.
 - + Independent weighing system

Each track is individually completed with high-precision weighing. Single-row or whole-column elimination can be realized.
 - + Strong anti-interference ability

The weighing machine independently supports the design, to ensure that the vibration of sorting and transmission, to ensure the weighing accuracy. Positive pressure protection, to ensure that the dust does not enter the scale body, to avoid affecting the weighing accuracy.
 - + Easy cleaning and maintenance.

The whole can be drawn out, the standard configuration of 15.6-inch split ultra-thin display screen + industrial control machine, simple and friendly interface.
 - + Automatic removal.

Single column can be eliminated, and the subsequent automatic bag filling mechanism should be selected.

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
This equipment is used with multiple lines of charter flight wiring. Complete the bag redundant storage, steering, bag filling, second count and other functions.

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Adhering to concept of “all for customers, for all customers” , and with the mission of “creating value for customers” , TRUKING focuses on building first-class R&D, design, sales, manufacturing, service teams and service management system to offer high-efficient and high-quality products and services for the customers.

◎ Project Planning Support

TRUKING has a series of professional medical and pharmaceutical design institutes, product research institutes, a group of senior scientific research and technical personnel, aseptic pharmaceutical process verification experts, professors, etc., to provide scientific and effective technical planning support for customers' projects.

◎ Process Assurance Service

TRUKING launched the overall solution for medical & pharmaceutical equipment realizing the full automation and aseptic production process, intelligent detection and sorting of the drug, as well as computer verification assurance so as to produce assured products for market.

◎ Service In The Park

As an equipment supplier, TRUKING not only integrates industrial regulations and standards into product design and manufacturing, but also introduces the industrial regulations and standards into Truking industrial park in the form of lectures and exchanges to interact and integrate with the technology and production of pharmaceutical enterprises.

◎ Equipment Management Consultation Service

TRUKING aims to satisfy the maximum benefited demands of the customer and to provide a package of equipment management services to realize minimization investment and maximization production output together with pharmaceutical factories.

◎ Adhere to the philosophy of "customer-centric", TRUKING creates the extraordinary service experience with the spirit of the craftsman, to create more value for customers.

◎ We are always ready to escort the normal operation of the customer's equipment.

- + Onsite equipment installation, commissioning, acceptance
- + Training for equipment process, operation and maintenance
- + Periodic door-to-door maintenance and return visit
- + Performance test and maintenance for machinery parts
- + Timely supply of spare parts
- + Product technology upgrading and secondary development
- + Equipment overhaul

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Truking Technology Industrial Park, Ningxiang National Economic and Technological Development Zone

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The company reserves the right to interpret the contents of this manual without prior notice. The pictures in this manual are for reference and subject to the actual conditions.

Apr. 2024 Edition