Blow-Fill-Seal Solution

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Equip global pharmaceutical manufacturing industry with intelligent wings.

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Founded in 2000, Truking Technology Limited has become a leading pharmaceutical equipment company in China and renowned worldwide. The company specializes in pharmaceutical equipment and integrated solutions and takes the initiative to push forward the R&D of smart pharmaceutical factories. It is a listed A-share company with multiple wholly-owned or holding subsidiaries, such as Germany-based Romaco Group, Truking Watertown, Sichuan Pharmaceuticals Design Institute, Truking Feiyun, Truking Ingenuity, Truking Micro-Sphere, Truking SUT, Truking Gene, Wachine, Truking Changxing, Truking Scientific Instrument, Truking Jingbang Engineering Technology, Truking Pitide Biotechnology, Truking Boyuan Intelligent Technology, Truking New Material, Truking Intelligent Robot etc. Truking has over 9000 employees around the globe, with a total asset of over 11 billion RMB.


Truking has two operational headquarters in Changsha, China, and Germany, with research and development institutions including Changsha Central Technology Research Institute, Europe Technology Research Institute, and Sichuan Pharmaceuticals Design Institute. It has established many platforms for technology innovation, such as the National Enterprise Technology Center, National Innovation Center, Post-Doctoral Scientific Research Station, and Academician Expert Workstation. By December 31, 2023, a total of 5,466 patent applications from China have been submitted, with 4,313 patents granted and 3,090 patents currently in force. Additionally, 54 PCT international patent applications were filed, resulting in 26 patents granted in multiple countries including the United States, Russia, India, South Korea, Germany, Indonesia, Japan, and Europe. The company has led the development of over 20 national industry product technical standards for this system. The products and services of the group have now reached over 180 countries and regions, with its international market share steadily increasing year by year.

吹 - 灌 - 封工艺概述 / Overview of blowing-filling-sealing process

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Truking BFS turnkey solution consists of six core functional modules: BFS equipment, liquid distribution system, sterilization and logistics system, lamp inspection and leak detection system, post packaging system and centralized control and scheduling system. Based on high efficiency and energy saving process equipment system, information system as the core, and automatic logistics as the link to achieve high efficiency, high quality BFS products, safe production.

工艺流程 / BFS Technical Process

挤塑 Extruding	<p>螺杆将粒子塑化连续挤出通过模头形成管坯。管坯内有无菌空气支撑。塑料热融和挤出的温度为 170~220 度，是经过实证的无菌挤出和去热源工艺。</p> <p>The Parison is melted and plasticized by extrusion screw and formed by extrusion head. The Parison is supported by sterile air. The temperature of plastic melted and extruded is 170~220 degrees</p>
成型 Blowing	<p>主模具合拢，通过真空和 / 或无菌压缩空气按照模具的形状，把塑料原料成型容器。</p> <p>The Parison is formed into container by vacuum and the sterile compressed air according to the mold shape when the main mold closing.</p>
灌装 Filling	<p>通过时间压力法将药液灌装至容器里，并排出无菌空气。药液经过的管路在灌装之前进行在线 CIP 和在线灭菌 SIP，确保药液安全性。</p> <p>Adopting Time-pressure methods to fill. Before filling, the line can be CIP & SIP to ensure the safety of liquid.</p>
封口 Sealing	<p>灌装结束后，将封口模具合闭，形成容器的顶部，并使瓶子密封。</p> <p>After filling, the sealing mold is closed to form the top of the container and the bottle is sealed.</p>
出瓶 Outfeeding	<p>瓶子封口结束后模具打开，已经完成成型、灌装并密封的瓶子被送出机器，进行切边。</p> <p>After sealing, the mold will be opened. The bottle will be transported to punching machine for separating the product and waste.</p>

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

Blow-fill-seal technology referred to as BFS. The technology is to complete the forming, filling and sealing of the product by reciprocating the die in the extrusion station and the filling station in a continuous operation process of one piece of equipment, all of which are completed under aseptic conditions, the whole cycle is completed within 12~14 seconds.

BFS technology has been used for about 40 years in the world, and it has been proved that the contamination rate can reach less than 0.1%, Based on the summary and analysis of the media filling data, the contamination rate of BFS system is 0.001%.

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+ Equipment capacity

The stable capacity of 24,000pcs/h (2ml), is the largest capacity of shuttle type BFS machine, with the mature technology and stable production.

+ Extruding system:

Designed from Truking (European) R&D Center, the extrusion die is manufactured with advanced machining, adopts the servo motor to drive, the Parison dynamic exhaust and fill gas to ensure the stability of the extruding system. Container wall thickness uniformity $\pm 0.05\text{mm}$.

+ Molding system:

A servo hydraulic system is used to control the PID closed-loop curve to reduce the assault and vibration of the equipment and avoid the wear of the mould. Servo swing cylinder drives high-speed mold shifting to reduce the risk of container exposure.

+ Filling system:

Adopting time-pressure methods to fill using high frequency solenoid valve control loading. Equipped with automatic identification technology of waste cylinders, the suspension uniformity deviation $\pm 0.5\%$, and the remaining oxygen $< 1\%$

+ Mold:

Adopting the precision manufacturing technology, manufactured by Truking mold center. and selects the imported non-toxic high-conductivity copper alloy for materials processing, to guarantee the service life of the mold and improve the product appearance quality.

+ Hydraulic systems:

Servo hydraulic pump station, PID adjustment of outlet pressure, closed-loop curve control of cylinder action. In the manual state, the hydraulic pump station runs at low speed, which is convenient for equipment debugging and maintenance and improves the safety protection of personnel.

+ Control system

Adopt SIEMENS motion controller and industrial computer to control, the action execution is more accurate implementation of electronic signature and electronic batch records to meet the data integrity verification requirements

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+ Equipment Capacity

The maximum stable production capacity of continuous BFS equipment designed by Truking (European) R&D Center is 33000pcs/h (2ml), and the utilization rate of plastic materials is up to 83%.

+ Extruding system

For the closed single extrusion parison, the die head is made of Europe-imported alloy materials that are processed by five-axis high-speed machine tool, adopting special coating for surface, and servo motor for extrusion screw. For Parison, adopting dynamic exhaust and air make-up so as to maintain the stability of extrusion system and uniformity of container wall thickness.

+ Molding System:

Multi-mold continuous equipment uses multiple pairs of molds installed on a continuous rotating chain, and the chain is made of alloy materials and thermally treated to guarantee its service life.

Single-mold continuous type equipment adopts the servo-hydraulic to open, close, lift and lowers the mold, which increases the closing force and improves the molding quality of the products.

+ Filling System:

The filling needles is installed inside the extrusion die head, protected by sterile air, and used to complete the filling of product in a closed parison. Time pressure method is used for filling, and high-frequency response electromagnetic valve island is adopted to control the filling volume, and the maximum filling accuracy can be up to $\pm 0.01\text{mL}$.

+ Forming Mold

The mold is manufactured by Truking mold center. The mold raw material is imported.

+ Control System

SIEMENS motion controller and industrial computer are used as control unit to realize electronic signature and electronic batch recording, so as to meet the requirements for data integrity verification.

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
Service Concept

400-9988-900

Adhering to concept of “all for customers, for all customers” , and with the mission of “creating value for customers” , TRUKING focuses on building first-class R&D, design, sales, manufacturing, service teams and service management system to offer high-efficient and high-quality products and services for the customers.

◎ Project Planning Support

TRUKING has a series of professional medical and pharmaceutical design institutes, product research institutes, a group of senior scientific research and technical personnel, aseptic pharmaceutical process verification experts, professors, etc., to provide scientific and effective technical planning support for customers' projects.

◎ Process Assurance Service

TRUKING launched the overall solution for medical & pharmaceutical equipment realizing the full automation and aseptic production process, intelligent detection and sorting of the drug, as well as computer verification assurance so as to produce assured products for market.

◎ Service In The Park

As an equipment supplier, TRUKING not only integrates industrial regulations and standards into product design and manufacturing, but also introduces the industrial regulations and standards into Truking industrial park in the form of lectures and exchanges to interact and integrate with the technology and production of pharmaceutical enterprises.

◎ Equipment Management Consultation Service

TRUKING aims to satisfy the maximum benefited demands of the customer and to provide a package of equipment management services to realize minimization investment and maximization production output together with pharmaceutical factories.

- ◎ Adhere to the philosophy of "customer-centric", TRUKING creates the extraordinary service experience with the spirit of the craftsman, to create more value for customers.
- ◎ We are always ready to escort the normal operation of the customer's equipment.

- + Onsite equipment installation, commissioning, acceptance
- + Training for equipment process, operation and maintenance
- + Periodic door-to-door maintenance and return visit
- + Performance test and maintenance for machinery parts
- + Timely supply of spare parts
- + Product technology upgrading and secondary development
- + Equipment overhaul

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