Integrated solution for separation and purification

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Equip global pharmaceutical manufacturing industry with intelligent wings

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Since 2012, the company has carried out various
types of chromatography, liquid preparation,
ultrafiltration and other equipment.
Development and manufacturing of equipment.

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Truking Ingenuity Biotechnology (Changsha) Co.,Ltd. (hereinafter referred to as Truking Ingenuity) is a holding subsidiary of Truking Technology Limited. Co-founded by Hunan United Ingenuity Technology Co., Ltd. (Stock code:Truking 300358)Truking Ingenuity positions itself as a comprehensive separation and purification solution provider for the global biopharmaceutical enterprises. It is located at the Truking Industrial Park in Ningxiang City, Hunan Province, China.

The core team of Truking Ingenuity originated from China Aerospace Science and Industry Corporation (CASIC), later was introduced in with international top talents from MILLIPORE and GE, and at the same time cooperated on technology and products with JNC Corporation, by jointly establish the China-Japan Technology Cooperation Process Lab. Truking Ingenuity was the earliest professional manufacturer of the separation and purification equipment of biological macromolecules in China. By now, Truking Ingenuity has enriched its product matrix to cover separation & purification, ultrafiltration, and Auto Inline Conditioning for labs from pilots to large-scale production. Process and auxiliary process development guidance is also available whenever the customer needs.

So far, Truking Ingenuity has established good cooperation with hundreds of companies home and abroad, with its products exported to many countries like Japan, India, Brazil, Russia, Cambodia, etc. The reliability, safety, stability, and convenience to use have been well recognized by our customers.

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

Truking Ingenuity Biotechnology (Changsha) Co.,Ltd. is committed to providing industry-leading chromatography equipment and technical services for global customers. At the same time, it has established a process development laboratory in Changsha in cooperation with JNC CORPORATION, to provide customers with overall solutions from chromatographic packing screening, purification process development, chromatographic equipment selection and customization.

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抗体纯化工艺/Actual process "process of Mabs"

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
The high specificity and affinity of antibody drugs that can bind to specific targets, and the easy humanization and fully humanization of the target protein sequence have led to the explosive growth of such biological products. This type of product has developed rapidly in the past 20 years, and more than 300 monoclonal antibodies are currently in clinical trials. Monoclonal antibodies have been approved for multiple indications such as malignant tumors, immunological diseases and rare diseases.

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The common antibody purification process is to first centrifuge and filter the harvested liquid from the reactor to remove cells and cell debris, and then pass the ProteinA chromatography to capture the antibody. After the captured antibody is virus inactivated by the low pH, it is passed through two ion exchange chromatography. Then filter the virus, and finally adjust the antibody concentration and change the liquid through ultrafiltration and diafiltration to obtain the antibody stock solution.

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Truking Ingenuity's new chromatographic process is based on the traditional purification process with one-step mixed mode chromatography instead of two-step ion exchange chromatography to achieve the separation and purification effect.

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Vaccines are microorganisms or their products that have been artificially treated. After being inoculated into the body, they can stimulate the immune system to produce antibodies and cells that specifically act on specific viruses or bacteria, so that the human body can develop resistance to the pathogen. Human virus vaccines mainly include live attenuated vaccines, inactivated vaccines, and viral vector vaccines. The production processes of these types of vaccines are similar.

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Virus separation and purification mainly uses the feature that virus particles are much larger than protein and peptide molecules, and is separated through membrane separation technology and chromatographic separation technology, and finally high-purity virus particles are obtained.

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The common virus particle purification process is to first remove the cell debris from the fermentation broth by microfiltration, then concentrate by Ultrafiltration to reduce the amount of sample and remove impurities such as HCP, and then pass Gel Filtration Chromatography and Anion Chromatography, After Ultrafiltration and Diafiltration concentration adjustment and liquid exchange, the stock solution is obtained.

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Truking Ingenuity's new purification process uses one-step affinity chromatography instead of gel filtration chromatography and anion chromatography on the basis of the traditional purification process to achieve the separation and purification effect of virus particles.

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Blood products are mainly a biologically active product prepared from healthy human plasma using biological processes or separation and purification technologies. Blood products include albumin, immunoglobulin, and various coagulation factors. In terms of production process, most blood products are purified by the traditional low-temperature ethanol method. In recent years, the production process of immunoglobulin and various coagulation factors has been gradually upgraded to use chromatography and ultrafiltration for separation and purification.

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Globulin is a serum protein that exists in human plasma and has an immune effect. It is also called immunoglobulin. The upgraded immunoglobulin production process is mostly used in ion chromatography, which mainly adopts two-step ion exchange column chromatography combining cation chromatography and anion chromatography, through two-step chromatography and ultrafiltration purification to meet the process requirements.

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Human coagulation factor VIII is a very important coagulation factor in human plasma. Human coagulation factor VIII absence or deficiency can cause varying degrees of hemophilia. Hemophilia A is a lifelong disease, and there is still no cure. Hemophilia A can only be treated by infusion of FVIII products. At present, the production process of human coagulation factor FVIII mainly adopts one-step ion exchange and ultrafiltration.

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Immunoglobulin purification process

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Coagulation factor VIII purification process

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Product introduction

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Macromolecule
purification series

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- + Adopt modular design for all valves, pumps and detectors, easy to dismantle and maintain;
- + Equip high accuracy binary piston pump with accurate flow flux and good repeatability;
- + Diversified configuration: with the most 8 channels buffer, 8 channels sample inlets, 8 channels outlets and 8 chromatography columns;
- + Optional automatic collector and sample loading pump to realize fully automatic sample loading and collecting;
- + Reliable computer managing system to ensure completeness, authenticity and traceability of the data;
- + Reasonable flow path design with small dead volume, stable liquid flow and accurate data measuring;
- + Powerful method editor to enable automatic sample loading, collection, monitoring, etc;
- + Multiple pipelines replaceable design with large range of flux and wide application;
- + Remote monitoring function to monitor the operation of equipment at any time anywhere.

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优势特点 / Performance advantages

- 结构紧凑, 堆叠式管路设计, 管路死体积小;
 - 专利阀反馈设计, 实时监测阀门开、关, 确保生产安全可靠;
 - 专利自主研发控制软件, 多代升级, 功能强大, 稳定可靠;
 - 简洁、直观的软件操作界面, 菜单化功能操作, 操作有序更易快速掌握;
 - 灵活而强大的方法编辑功能, 可基于时间、体积、柱体积进行自动流程生产;
 - 远程监控功能, 随时随地监测设备运行状态。
- + With compact structure, stack type design and small dead volume;
 - + With patented valve feedback design, which could realize real-time monitoring for valve opening and closing, to ensure the safety production;
 - + Stable, reliable and powerful functional self independent research and development patented control software with multiple generation upgrade;
 - + Concise and intuitive software operation interface with menu function operation to enable orderly operation and easy to fast master;
 - + Flexible and powerful method editing function enables automatic process production based on time, volume and column volume;
 - + Remote monitoring function to monitor the running status of equipment at any time anywhere.

产品型号 Product model	AP400S(D)	AP800S(D)	AP1600S(D)	AP2500S(D)
流量(L/min) Flow rate (L/min)	2.5~40	5~80	10~160	20~250
主管路尺寸 Main piping size	DN25(1")	DN40(1.5")	DN50(2")	DN65(2.5")
标配气泡陷阱体积(ml) Standard Bubble Trap Volume(ml)	4000	8000	16000	25000
流量精度 Flow Accuracy	检测精度±0.5%, 系统流量精度±2% Detection accuracy is ±0.5%, flow accuracy is ±2%			
压力传感器范围 Pressure Sensor Range	测量范围:0~10 bar, 精度 ≤ 0.5% Measuring range: 0~10 bar, accuracy ≤ 0.5%			
pH 范围 pH Range	测量范围:0~14 ;最小分辨率0.01;精度±0.06 Measuring Range: 0~14; min. resolution: 0.01; accuracy ±0.06			
电导范围 Conductance Range	测量范围Measuring Range:1 μS/cm - 300 mS/cm; 1 μS/cm to 100 mS/cm, 精度Accuracy±3% 100 to 300 mS/cm, 精度Accuracy±5%			
UV范围 UV Range	波长范围:200~400 nm或190~700 nm, 检测范围:0~6AU, 精度±1nm; 可选:固定和可调, 单通道, 双通道, 四通道 Wavelength Range:200~400 nm/190~700 nm, Measuring Range: 0~6AU, accuracy ±1nm; Optional: FIX or DAD, 1 channel, double channel, four channel			
防护级别 Protection class	IP54			
接液材质 Wetted material	316L、PP、EPDM、高硼硅玻璃 high borosilicate glass、PEEK			
工作温度范围 Operating Temperature Range	2~40℃			
最大操作压力 Max. Operating Pressure	6 bar			
供电 Power Supply	AC, 200~240V, 50/60Hz		3AC, 340V~420V, 50/60Hz	
功率 Power(KW)	3 (4)	4 (7)	5 (8)	6 (11)
尺寸(L×W×H) Dimensions (L×W×H)	1300×740×1200 mm 琴式 1500×700×1770 mm 旋臂式	2340×1350×2300 mm 旋臂式Swivel	2475×1580×2650 mm 旋臂式Swivel	2475×1580×2650 mm 旋臂式Swivel
设备重量(kg) Weight (kg)	550/900	1200	1600	2000

层析系统 PURIFICATION

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MULTIPLE COLUMN
CONTINUOUS PURIFICATION

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- + Through the series connection of chromatography columns, full load loading is realized, and the utilization rate of fillers is improved;
- + Reduce the liquid preparation and the corresponding storage tank design, save 50% of the buffer consumption, and miniaturize the equipment;
- + Improve the purity and yield of the product.

主要技术参数 / Main technical parameters

产品型号 Product model	MCC010T	MCC030T	MCC100T	MCC200T
流量精度 Flow Accuracy	0.02~1	0.05~3	0.3~10	0.5~20
主管路尺寸 Main piping size	DN08(1/4")	DN10(3/8")	DN15(1/2")	DN20(3/4")
标配气泡陷阱体积(mL) Standard Bubble Trap Volume(ml)	105	300	1500	1500
出入口数量 Number of Inlets and Outlets	每组1~8路 1~8 Channels per group			
流量精度 Flow Accuracy	检测精度±0.5%, 系统流量精度±2% Detection accuracy is ±0.5%, flow accuracy is ±2%			
压力传感器范围(bar) Pressure Sensor Range(bar)	测量范围:0~10 bar, 精度 ≤ 0.5% Measuring range: 0~10 bar, accuracy ≤ 0.5%			
pH 范围 pH Range	测量范围:0~14 ;最小分辨率0.01;精度±0.06 Measuring Range: 0~14; min. resolution: 0.01; accuracy ±0.06			
电导范围 (ms/cm) Conductance Range (ms/cm)	测量范围Measuring Range:1 μS/cm - 300 mS/cm; 1 μS/cm to 100 mS/cm, 精度Accuracy±3% 100 to 300 mS/cm, 精度Accuracy±5%			
UV范围 UV Range	波长范围:200~400 nm或190~700 nm, 检测范围:0~6AU, 精度±1nm; 可选: 固定和可调, 单通道, 双通道, 四通道 Wavelength Range:200~400 nm/190~700 nm, Measuring Range: 0~6AU, accuracy ±1nm; Optional: FIX or DAD, 1 channel, double channel, four channel			
密封及保护级别 Protection class	IP54			
接液材质 Wetted material	316L, PP, EPDM、高硼硅玻璃 high borosilicate glass、PEEK			
最大操作压力 Operating pressure max.	2~40°C			
工作温度范围 Operating Temperature Range	5 bar			
供电 Power Supply	AC, 200~240V, 50/60Hz			
功率 Power(KW)	1	1.5	3	4
尺寸(L×W×H) Dimensions (L×W×H)	530×520×740 mm		1230×790×1920 mm	
设备重量(kg) Weight (kg)	260	265	290	300

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- + The series design of products enriches the choices of customers and can be customized according to customer requirements;
- + Graphical interface display, full touch operation;
- + Humanized electrical components and circuit layout;
- + Device does not require CIP/SIP;
- + Pipe fittings are easy to disassemble;
- + Self-developed control software, in line with FDA and cGMP;
- + Newly developed with GE Unicorn software as a template for flexible method editing.

主要技术参数 / Main technical parameters

产品型号 Product model	APU010S(D)	APU030S(D)	APU100S(D)	APU200S(D)	APU400S(D)
流量(L/min) Flow rate (L/min)	0.02~1	0.05~3	0.3~10	0.5~20	2.5~40
主管路尺寸 Main piping size	6.4 mm	9.5 mm	12.7 mm	15.9 mm	25.4 mm
系统泵 System Pump	耐压 4bar; 泵头Gamma射线消毒, 一次性安装; 接液材质: PP、TPE、EPDM Resistant to pressure 4 bar, pump head gamma radiation sterilization, one-time installation; wetted materials: PP, TPE, EPDM				
一次性更换器件 Disposable Size Parts	泵头、管道、pH 探头、电导率探头、压力探头、紫外流通池、气泡陷阱、过滤器 (选配) Pump head, pipeline, pH detector, conductance detector, pressure detector, UV flow cell, bubble trap, filter (optional)				
流量精度 Flow Accuracy	±2%				
压力传感器范围 (bar) Pressure Sensor Range(bar)	-0.48~5.2 bar				
pH 范围 pH Range	测量范围: 3~10; 最小分辨率: 0.01; 精度 ±0.06 Measuring Range: 3~10; min. resolution: 0.01; accuracy ±0.06				
电导传感器 Conductance Sensor	测量范围 Measuring Range: 1 μS/cm - 300 mS/cm; 1 μS/cm to 100 mS/cm, 精度 Accuracy ±3% 100 to 300 mS/cm, 精度 Accuracy ±5%				
UV 范围 UV Range	波长范围: 200~400 nm 或 190~700 nm, 检测范围: 0~6AU, 精度 ±1nm; 可选: 固定和可调, 单通道, 双通道, 四通道 Wavelength Range: 200~400 nm/190~700 nm, Measuring Range: 0~6AU, accuracy ±1nm; Optional: FIX or DAD, 1 channel, double channel, four channel				
密封及保护级别 Protection class	IP54				
接液材质 Wetted material	316L、PP、EPDM、高硼硅玻璃 high borosilicate glass				
工作温度范围 Operating Temperature Range	2~40°C				
最大操作压力 Operating pressure max.	4bar (波纹管路), 1 bar (硅胶管路) 4 bar (corrugated pipeline), 1 bar (silicone pipeline)				
供电 Power Supply	AC, 200~240V, 50/60Hz				
功率 Power(KW)	0.8	1	1.5	2	3
尺寸 (L × W × H) Dimensions (L × W × H)	1850 × 1100 × 2100 mm				
设备重量 (kg) Weight (kg)	240	240	250	260	400

层析系统 PURIFICATION

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AUTO INLINE CONDITIONING
CHROMATOGRAPHY

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- + Only high concentration stock solution is needed, which saves a lot of floor space and space;
- + Precise configuration of Buffer;
- + High repeatability of different batches;
- + Automatic control by computer, built-in 6-way dynamic controller, the configuration of different batches can be perfectly reproduced.

主要技术参数 / Main technical parameters

产品型号 Product model	ICC030Q	ICC100Q	ICC200Q	ICC400Q	ICC800Q	ICC1600Q	ICC2500Q
流量(L/min) Flow rate (L/min)	0.05~3	0.3~10	0.5~20	2.5~40	5~80	10~160	20~250
主管路尺寸 Main piping size	DN10(3/8")	DN15(1/2")	DN20(3/4")	DN25(1")	DN40(1.5")	DN50(2")	DN65(2.5")
标配气泡陷阱 体积(ml) Standard Bubble Trap Volume(ml)	300	1500	1500	4000	8000	1600	25000
流量精度 Flow Accuracy	检测精度±0.5%, 系统流量精度±2% Detection accuracy is ±0.5%, flow accuracy is ±2%						
压力传感器范围 Pressure Sensor Range	测量范围:0~10 bar, 精度 ≤ 0.5% Measuring range:0~10 bar, accuracy ≤ 0.5%						
pH范围 pH Range	测量范围:0~14; 最小分辨率0.01; 精度±0.06 Measuring Range: 0~14; min. resolution: 0.01; accuracy ±0.06						
电导范围 Conductance Range	测量范围Measuring Range:1 μS/cm - 300 mS/cm; 1 μS/cm to 100 mS/cm, 精度Accuracy±3% 100 to 300 mS/cm, 精度Accuracy±5%						
UV范围 UV Range	波长范围:200~400 nm或190~700 nm, 检测范围:0~6AU, 精度±1nm; 可选: 固定和可调, 单通道, 双通道, 四通道 Wavelength Range:200~400 nm/190~700 nm, Measuring Range: 0~6AU, accuracy ±1nm; Optional: FIX or DAD, 1 channel, double channel, four channel						
防护级别 Protection class	IP54						
接液材质 Wetted material	316L、PP、EPDM、高硼硅玻璃 High borosilicate glass、PEEK						
工作温度范围 Operating Temperature Range	2~40°C						
最大操作压力 Operating pressure max.	6 bar						
供电 Power Supply	AC, 200~240V, 50/60Hz			3AC, 340V~420V, 50/60Hz			
功率 Power(KW)	1.5	2	4	5	10.5	12	16
尺寸(L×W×H) Dimensions (L×W×H)	1900×930×1650 mm			1600×800 ×1850 mm	2340×1350 ×2300 mm	2475×1580×2650 mm	
设备重量(kg) Weight (kg)	550	705	725	950	1200	1600	2000

层析系统 PURIFICATION

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SMART PURIFICATION

优势特点 / Performance advantages

- ◎ 结构紧凑简洁,小巧轻便,运转方便;
 - ◎ 管路多种搭配,灵活选择,拆卸、维护方便;
 - ◎ 搭配手动四通阀,长寿命,低温不漏,高温可SIP;
 - ◎ 手动四通阀可实现气泡陷阱旁路和在线,同时能实现层析柱正、反向和旁路;
 - ◎ 拥有调速模式、流量模式、压力模式三种泵控模式,实现多样化操作;
 - ◎ 自主研发控制软件,简洁、直观的软件操作界面,菜单化功能操作,操作有序更易快速掌握;
 - ◎ 集数据采集、显示、存储于一身,符合GMP要求。
- + Compact and simple structure design, small in size and convenient for operation;
 - + Multiple choices for pipelines, convenient for dismantling and maintenance;
 - + Equip patented manual four-way valve with long life span, no leakage under low temperature and SIP is available in high temperature;
 - + Manual four-way valve to realize the bypass and on-line of bubble trap, as well as the forward, reverse and bypass of chromatographic column;
 - + Three pump control modes: speed regulation mode, flow mode and pressure mode are available to realize diversified operation;
 - + Self research and developed control software with concise and intuitive software operation interface and menu operation to enable orderly operation and easy to master quickly;
 - + Integrated data collection, display, storage as a whole meeting the requirements of GMP.

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- + Adopting ultrafiltration membrane clamp matching the well-known brand ultrafiltration membrane standards;
- + Product design is in series and modular, which can be customized according to customer process requirements;
- + All structures adopt sanitary structure design and pipeline is stacked type with features of compact structure and the smallest dead volume of pipeline in the industry;
- + Configured with touch all-in-one machine, available for graphical interface display and full touch operation;
- + Self developed control system which is comply with the requirements of FDA and cGMP;
- + The software interface is concise and intuitive with menu function setting to enable more orderly operation;
- + Reliable computer managing system to ensure completeness, authenticity and traceability of the data;
- + Remote monitoring function to monitor machine running status at any time anywhere.

主要技术参数 / Main technical parameters

产品型号 Product model	AF0050	AF0100	AF0200	AF0300	AF0400	AF0800	AF1000
流量(L/min) Flow rate (L/min)	5~50	10~100	20~200	30~300	40~400	80~800	100~1000
主管路尺寸 Main piping size	DN25(1")	DN40(1.5")	DN50(2")	DN65(2.5")	DN80(3")	DN100(4")	DN100(4")
最大膜面积(m ²) Max. Membrane (m ²)	5	10	20	30	40	80	100
压力传感器范围 Pressure Sensor Range	测量范围:0~10 bar, 精度 ≤ 0.5% Measuring range: 0~10 bar, accuracy ≤ 0.5%						
pH范围 pH Range	测量范围:0~14 ;最小分辨率0.01;精度±0.06 Measuring Range: 0~14; min. resolution: 0.01; accuracy ±0.06						
电导范围 Conductance Range	测量范围Measuring Range:1 μS/cm - 300 mS/cm; 1 μS/cm to 100 mS/cm, 精度Accuracy±3% 100 to 300 mS/cm, 精度Accuracy ±5%						
UV范围 UV Range	波长范围:200~400 nm或190~700 nm, 检测范围:0~6AU, 精度±1nm; 可选: 固定和可调, 单通道, 双通道, 四通道 Wavelength Range:200~400 nm/190~700 nm, Measuring Range: 0~6AU, accuracy ±1nm; Optional: FIX or DAD, 1 channel, double channel, four channel						
控制模式 Control Model	恒TMP、恒压、恒流、恒速、恒重自动补液、自动水通量测试 Constant TMP, constant pressure, constant current, constant speed, constant weight automatically replenishes the liquid, ultrafiltration Membrane water flux testing						
防护级别 Protection class	IP54						
接液材质 Wetted material	316L、PP、EPDM、高硼硅玻璃 High borosilicate glass、PEEK						
工作温度范围 Operating Temperature Range	2~40°C						
最大操作压力 Operating pressure max.	6 bar						
供电 Power Supply	3AC, 340V~420V, 50/60Hz						
功率 Power(KW)	3	6	6	9	9	9	17
尺寸(L × W × H) Dimensions (L × W × H)	1800 × 800 × 1800 mm	1500 × 1100 × 1900 mm	2100 × 1500 × 2300 mm	2500 × 1500 × 2500 mm	2500 × 1500 × 2500 mm	2000 × 2000 × 3500 mm	2200 × 2200 × 3500 mm
设备重量(kg) Weight (kg)	700	950	1500	1700	1800	2300	2500

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EASY TFF

优势特点 / Performance advantages

- ⊙ 设备整体设计简洁紧凑，方便维护和转移；
 - ⊙ 配置手动单阀，可手动控制调节TMP；
 - ⊙ 可匹配多种品牌膜包的夹具；
 - ⊙ 路最短化设计，系统无残留，清洗更简单；
 - ⊙ 全部卫生型结构设计，符合FDA要求的材质。
- + The overall design of the equipment is simple and compact, which is convenient for maintenance and transfer;
 - + Equipped with manual single valve, which can manually control and adjust TMP;
 - + Holder can match various brands of Membranes; Shortest pipeline design, no residue in the system, easier cleaning;
 - + All hygienic structural design, in line with FDA requirements.

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- + The series design of products enriches the choices of customers and can be customized according to customer requirements;
- + The control cabinet is made of 304/316L stainless steel, which is waterproof, rust-proof and easy to clean;
- + Patented valve feedback design, real-time monitoring of valve opening and closing to ensure safe production;
- + Humanized and beautiful electrical components and circuit layout;
- + The software interface is simple and intuitive, and the menu-based function settings make the operation more orderly;
- + Stable and reliable data management, real-time recording and traceability, no data loss after power failure;
- + All hygienic structural design, in line with FDA requirements.


主要技术参数 / Main technical parameters

产品型号 Product model	DF0020	DF0040	DF0080	DF0160
流量(L/min) Flow rate (L/min)	1~20	5~40	10~80	20~160
最大膜面积(m ²) Max.Membrane (m ²)	两级过滤, 每级5m ² Two-stage filtration, each stage is 5m ²	两级过滤, 每级10m ² Two-stage filtration, each stage is 10m ²	两级过滤, 每级20m ² Two-stage filtration, each stage is 20m ²	两级过滤, 每级40m ² Two-stage filtration, each stage is 40m ²
适配过滤膜品牌 Membrane brand	Millipore、Sartorius、3M、PALL			
压力传感器范围 Pressure Sensor Range	测量范围Measuring range:0~10 bar, 精度accuracy ≤ 0.5%			
电导范围 Conductance Range	测量范围Measuring Range:1 μS/cm - 300 mS/cm; 1 μS/cm to 100 mS/cm, 精度Accuracy ± 3% 100 to 300 mS/cm, 精度Accuracy ± 5%			
选配: 光学浊度计 Optional: optical turbidimeter	双通道Dual channel, 0~100NTU			
控制模式 Control Model	恒压、恒流、恒速 Constant pressure, constant current, constant speed			
防护级别 Protection class	IP54			
接液材质 Wetted material	316L、PP、EPDM、高硼硅玻璃 High borosilicate glass、PEEK			
工作温度范围 Operating Temperature Range	2~40°C			
最大操作压力 Operating pressure max.	6 bar			
供电 Power Supply	AC, 200~240V, 50/60Hz		3AC, 340V~420V, 50/60Hz	
功率 Power(kw)	2	3	4	5
尺寸(L×W×H) Dimensions (L×W×H)	1400×600×1600 mm		1600×800×1750 mm	
设备重量(kg) Weight (kg)	300	400	550	600

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- + The series design of products enriches the choices of customers and can be customized according to customer requirements;
- + Graphical interface display, full touch operation;
- + Humanized electrical components and circuit layout;
- + Device does not require CIP/SIP;
- + It can be equipped with integrity test function;
- + Automatically test water flux;
- + Compatible with various brands of membrane;
- + Pipe fittings are easy to disassemble;
- + Self-developed control software, in line with FDA and cGMP;
- + New research and development using GE Unicorn software as a template to achieve flexible method editing.

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- + Support customized design, accurate and high efficiency;
- + Using computer instead of human to calculate and adjust solution preparation;
- + Available for three channels, four channels or multiple channels online solution preparation requirements according to customer requests;
- + Real time online detection to enable accurate parameter feedback and control, which can reduce the labor cost and manual operating time and error to improve efficiency;
- + To reduce enterprise costs in great extent and maximize convenience;
- + Patented valve feedback design monitoring the open and close of the valve in real time to ensure safe production;
- + Optimal PID control, scientific logic design;
- + Safe and reliable computer managing system to ensure completeness, authenticity and traceability of the data.

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- + Dual motor drive design, more reliable compared with single motor drive, available for synchronous monitoring and comparison to avoid motor failure.
- + Patented structural mechanics design to simplify equipment structure and strengthen whole equipment strength, enabling the equipment with less weight, high strength, small volume and lower height.
- + With motor concealed design and higher hygienic grade full PTFE dust cover design, the machine is more reliable and having a good appearance.
- + Patented fluid simulation distributor design, with features of uniform distribution, high column efficiency, good symmetry and reproducibility.
- + Patented innovative sieve plate fixed structure design, convenient for air exhaust and changeover with very low cost.
- + Patented column tube dismantling structure design, which greatly improves cleaning and maintenance efficiency and small space is needed.
- + Patented positive pressure column packing method, the column bed height designed to be higher to prevent the generation of slurry bubbles and greatly improve the utilization rate of filler. The filler consumption is the lowest in the industry.

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Note: The material of the column tube is divided into two types: acrylic and stainless steel. Acrylic is divided into two types: imported from the UK and imported from Japan, and relevant material certificates can be provided.

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- + Adopting patented column head sealing method which is more reliable than pneumatic sealing and lasting longer than spring-driven sealing;
- + Patented column head side-sliding and deflection design, convenient for dismantling and maintenance;
- + Equipped with a patented anti-low temperature leakage four-way valve, no fear of temperature difference, safe and reliable;
- + 316L stainless steel mesh design, with features of higher strength and no deformation;
- + Minimum sealing structure, with features of reliable performance and long life span;
- + Reliable sieve plate fixing method to maximize the effective distribution area;
- + Patented fluid simulation distributor design to ensure optimal distribution effect;
- + Multiple way fixed radial and axial directions, professional axial seal design, safe and reliable;
- + Complete verification documents and services meet GMP requirements.

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Description: The manual chromatography column can be equipped with a column extender to increase the height of the column, and can be customized according to customer needs.

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- + Automatic guided operation mode: there are guided instructions for column packing, unpacking, cleaning, and disassembly, which simplifies the operation and makes it easy to get started;
- + Flexible selection of various operation modes: the column can be operated with absolute positioning and relative movement, and at the same time, the manual flow path control can be performed in the automatic guided operation mode, which is simple and convenient;
- + Various stability and safety protection monitoring functions: The packing station has various monitoring settings, and the corresponding information can be viewed when the equipment alarms, so as to fully protect the safe and stable operation of the equipment.

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- + Configured with control and display system, it can realize the functions of packing slurry, valve automatic switching, slurry preservation solution system replacement and real-time slurry weight display;
- + It can be used with the automatic packing station to realize the slurry, automatic packing and unpacking of the full-automatic chromatographic column media;
- + Double screw impeller design with features of low rotary speed, low shear and fast mixing;
- + Depth air inlet design to pre-mix the media through clean compressed air;
- + Equip automatic weighing system to display real time weight inside of the tank accurately;
- + Double removable mesh, easy for alcohol replacement and mesh cleaning and maintenance;
- + Double spray ball design, online CIP cleaning, no dead corner.

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Description: Customized design can be made according to customer needs.

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Harmonic delivery pumps can be used in many manufacturing fields such as laboratory research and development, pharmaceutical industry, food industry, chemical industry, etc. They can gently, safely and reliably deliver shear-sensitive aqueous solutions and biological products.

主要技术参数 / Main technical parameters

产品型号 Product model	UFP018-S	UFP180-S	UFP1000-S
最大流量 (L/min) Maximum flow (L/min)	0.3	3	16
最小流量 (L/min) Minimum flow (L/min)	0.005	0.05	0.25
压力范围(bar) Pressure range (bar)	0-2	0-2	0-2
最大转速(rpm) Maximum speed (rpm)	1800	1000	1000
工作温度(°C) Working temperature(°C)	2-40	2-40	2-40
连接器 Connector	M8 (内螺纹) 或宝塔接头 M8 internal thread or pagoda joint	Tc卡箍Clamp (φ25)	TC卡箍Clamp (φ25)
接液材质 Wetted material	卫生级PP, PEEK, BPT管 (品牌: 圣戈班, 符合FDA认证) Sanitary PP, PEEK, BPT pipes (Brand: Saint Gobain, FDA compliant)	卫生级PP, BPT管 (品牌: 科尔帕默, 符合FDA认证) Sanitary PP, BPT pipe (brand: Cole-Parmer, FDA compliant)	卫生级PP, BPT管(品牌: 圣戈班, 符合FDA认证), 密封圈 (EPDM70) Sanitary PP, BPT pipe (brand: Saint-Gobain, FDA compliant), sealing ring (EPDM70)
泵腔体积(ml) Pump chamber volume (ml)	0.6	6.6	35
功率(W) Power (W)	100	100	1000
电源(V) Power (V)	AC220V	AC220V	AC220V
重量(Kg) Weight(Kg)	泵头Pump head: 0.6 整机Whole machine: 4	泵头Pump head: 1.4 整机Whole machine: 4.8	泵头Pump head: 8.1 整机Whole machine: 29
其他说明 Other explanations	CIP、SIP	CIP、SIP	CIP、SIP
长宽高 Length width height	泵头Pump head 71*50*124 整机Whole machine 222*181*204	泵头Pump head 129*93*182 整机Whole machine 236*181*204	泵头Pump head 285*130*380 整机Whole machine 305*328*448

优势特点 / Performance advantages

- | | |
|-------------------|---|
| ◎ 流态均匀稳定, 剪切力低; | + The flow status is even and stable with low shear force; |
| ◎ 流量精度高, 可正反输送; | + High flow flux accuracy, available to transmitted forwardly and backwardly; |
| ◎ 寿命超长, 免维护; | + Long life span, maintenance free; |
| ◎ 安全性高, 不会超压损坏; | + High safety, no damage under over voltage situation; |
| ◎ 自关断, 防虹吸, 防自流; | + Self closing, anti siphon and anti gravity flow. |
| ◎ 卫生等级高, 符合GMP标准。 | + High hygienic level, comply with GMP standard; |

产品型号 Product model	UFP180SU-S	UFP1000SU-S	UFP180M	UFP1000M
最大流量 (L/min) Maximum flow (L/min)	3	16	3	16
最小流量 (L/min) Minimum flow (L/min)	0.05	0.5	0.05	0.25
压力范围(bar) Pressure range (bar)	建议不带压 No pressure is recommended		0-6	0-6
最大转速(rpm) Maximum speed (rpm)	1000	1000	3000	2000
工作温度(°C) Working temperature(°C)	2-40	2-40	2-40	2-40
连接器 Connector	Y型三通Y type tee (Y064)	Y型三通Y type tee (Y127)	Tc卡箍TC clamp (φ25)	
接液材质 Wetted material	卫生级软管 Sanitary hose		卫生级PP, PEEK, 密封圈 (EPDM70) BPT管(品牌: 圣戈班, 符合FDA认证) Sanitary PP, PEEK, sealing ring (EPDM70) BPT pipe (brand: Saint-Gobain, FDA-compliant)	
泵腔体积(ml) Pump chamber volume (ml)	6.6	30	6.9	43
功率(W) Power (W)	100	1000	200	1000
电源(V) Power (V)	AC220V	AC220V	AC220V	AC220V
重量(Kg) Weight(Kg)	泵头Pump head: 3.4 整机Whole machine: 7	泵头Pump head: 14.1 整机Whole machine: 27	/	/
其他说明 Other explanations	易换管 Easy change tube		CIP, SIP	CIP, SIP
长宽高 Length width height	泵头Pump head 126*99*217	泵头Pump head 205*210*335	/	/
	整机Whole machine 245*193*275	整机Whole machine 360*300*385	/	/

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Remark: Customized design is available according to customer requests. For more information, please contact Truking Ingenuity Biotechnology (Changsha) Co., Ltd.

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+ Patented structure design to prevent low temperature leakage;
+ Precision grinding of sealing surface, surface finish within $ra0.02$, reliable sealing.

产品型号 Product model	MV020404L	MV020406L	MV020408L	MV020410L
通径 (mm) Drift diameter (mm)	4.57	6	7.75	9.4

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+ Non rotating manual homogenization, flange hole structure, no damage to media, high homogenization efficiency;
+ The material of PP meets the requirements of FDA.

产品型号 Product model	M001	M002
直径 (mm) Diameter (mm)	80	140

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+ Suitable for millipore, sartorius, Pall membrane;
 + 316L material, small residual volume, no cleaning dead corner;
 + Precision grinding of surface and inner hole, roughness Ra0.4;
 + The surface treatment technology of combining electropolishing and passivation in Germany technology.

产品型号 Product model	CH003	CH025	CH050	CH200
最大可夹持膜包面积 (m ²) Maximum clampable Membrane area(m ²)	0.3	2.5	5	20
进料 (inch) Feed (inch)	1/2"	3/4"	1"	1.5"
回流 (inch) Reflow (inch)	1/2"	3/4"	1"	1.5"
透过 (inch) Penetration (inch)	3/8"	1/2"	3/4"	1"
重量 (kg) Weight (kg)	6	18	32	50
长*宽*高(mm) Length*width*height(mm)	210*135*280	260*250*200	390*300*220	830*335*248
特点 Characteristics	表面光洁度Ra 0.4 电抛；接液材质316L；过程连接采用Tri-Clamp 连接；Millipore、Sartorius、Pall 盒式平板膜 Surface finish Ra 0.4 electropolished; wetted material 316L; process connection using Tri-Clamp Ligation; Millipore, Sartorius, Pall cassette plate membranes			

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All the chromatography system, ultrafiltration system and preparation system of Truking Ingenuity adopt self developed control software - PuriCAT software from laboratory level to production level with good continuity.

Software with multiple generation upgrade, easy for operation, convenient for usage, good stability and easy to master.

Available to install on any personal computer(PC) which has windows operation system, compatible with keyboard or touch screen.

PuriCAT is the software trade mark of Truking Ingenuity Biotechnology (Changsha) Co., Ltd.

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In addition to PuriCAT software that supports stand-alone operation, all purification equipment supports DCS batch automation solution of the whole plant.

The PLC/DCS programs are written in accordance with ISA88 standard to decouple control programming and process programming, so as to realize batch control and full automation of production process.

In addition to PuriCAT software that supports stand-alone operation, all At present, there are Emerson DeltaV, Rockwell PlantPAX, SIEMENS PCS7 and other implementation cases, and can be connected to any plant DCS system, MES system or SCADA system through network interface or OPC.

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
The main interface of Puricat software has divided into 3 areas: System function area, process control area, status display area;

The whole interface is concise and intuitive with menu operation mode, which is easy to master;


The process control area is responsible to control the valve, pump and detector directly and to check all sensors value and status;

Different system control area with different flow path model, the other interfaces remain the same.


生产级层析系统软件界面 | Production Chromatography System Software Interface

生产级超滤系统软件界面 | Production Ultrafiltration System Software Interface

在线配液层析系统软件界面 | Automatic Inline Conditioning Chromatography System Software Interface

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
- 拥有泵速、流量、压力、TMP、pH/电导等多种控制模式。
- Own multiple control modes, like control of pump speed, flow rate, pressure, TMP, pH and conductance, etc.

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

- 虚实标线, 查看详细信息; 曲线对比; 自动计算柱效等。
- Fictitious and actual marking. Dotted and solid lines to check the detailed information, curve comparison, automatic calculation of column efficiency, etc.

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
Truking Micro-Sphere Biotechnology (Changsha) Co., Ltd. (referred to as Truking Micro-Sphere) is a subsidiary of Truking Science and Technology Holdings, established in September 2021. The company focuses on the R&D and production of high-precision, high-performance and high-value-added chromatography microspheres, and provides international-quality products and services for biopharmaceutical companies and R&D institutions such as antibodies, vaccines, blood products, gene/cell therapy, and recombinant proteins.

Truking Microsphere Company can precisely manufacture a variety of chromatography media, including agarose, dextran and silica products. Truking Microsphere Company specializes in the functional modification of the surface of microspheres, the coupling of ligands and scaffolds, and has large-scale production capacity. Products are mainly used in the field of separation and purification of natural organics, antibodies, vaccines, blood products, gene/cell therapy, recombinant proteins and other biological products.

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- + The materials have biosafety and are widely used in the purification of biomolecules in food, cosmetics, medicine and other fields;
- + The product has excellent mechanical properties and can withstand pressure and alkali treatment;
- + There are many kinds of products, which can meet the needs of different customers such as laboratory process development and factory scale production.

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We have deep expertise, regulatory knowledge and many years Practical experience to develop a verification strategy for you in depth.

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
Complete verification document system, based on customer-specific verification Certificate documents (Chinese/English), in line with EU, US, ChinaGMP, Australian TGA, WHO and other global GMP regulations Require.

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0731-85051996 

Service Concept

Adhering to concept of "all for customers, for all customers", and with the mission of "creating value for customers", TRUKING focuses on building first-class R&D, design, sales, manufacturing, service teams and service management system to offer high-efficient and high-quality products and services for the customers.

Project Planning Support

TRUKING has a series of professional medical and pharmaceutical design institutes, product research institutes, a group of senior scientific research and technical personnel, aseptic pharmaceutical process verification experts, professors, etc. to provide scientific and effective technical planning support for customers' projects.

Process Assurance Service

TRUKING launched the overall solution for medical & pharmaceutical equipment realizing the full automation and aseptic production process, intelligent detection and sorting of the drug, as well as computer verification assurance so as to produce assured products for market.

Service In The Park

As an equipment supplier, TRUKING not only integrates industrial regulations and standards into product design and manufacturing, but also introduces the industrial regulations and standards into Truking industrial park in the form of lectures and exchanges to interact and integrate with the technology and production of pharmaceutical enterprises.

Equipment Management Consultation Service

TRUKING aims to satisfy the maximum benefited demands of the customer and to provide a package of equipment management services to realize minimization investment and maximization production output together with pharmaceutical factories.

■ Adhere to the philosophy of "customer-centric", TRUKING creates the extraordinary service experience with the spirit of the craftsman, to create more value for customers.

■ We are always ready to escort the normal operation of the customer's equipment.

- + Onsite equipment installation, commissioning, acceptance
- + Training for equipment process, operation and maintenance
- + Periodic door-to-door maintenance and return visit
- + Performance test and maintenance for machinery parts
- + Timely supply of spare parts
- + Product technology upgrading and secondary development
- + Equipment overhaul

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Add:Truking Technology Industrial Park,Ningxiang National
Economic and Technological Development Zone

产品技术改进恕不另行通知, 本公司保留对本手册内容的最终解释权, 本册图片仅供参考, 以实物为准。2023年3月版
The company reserves the right to interpret the contents of this manual without prior notice. The pictures in this manual are for reference and subject to the actual conditions.
March 2023 Edition